
Gravity 6 Long Term Care Pressure Redistribution Mattress

USER MANUAL

Item # 15670

USER: BEFORE using this product read this manual, containing important information on cleaning and maintenance, and save for future reference.

1 General

1.1 Symbols

Signal words are used in this manual and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING

Warning indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

Gives useful tips, recommendations and information for efficient, trouble-free use.

2 Cleaning and Disinfection Instructions

WARNING

Mattresses **MUST** be cleaned and disinfected regularly after each patient use. Because of the potential risk of infectious exposure, cleaning with the patient on the bed is not recommended. All equipment should be visually inspected. Any item that is visibly soiled with the patient's blood or other bodily fluids should be properly cleaned or removed. Staff members should treat all soiled bedding using universal precaution: as if it were contaminated with pathogenic microorganisms. Staff members should wear appropriate Personal Protective Equipment when cleaning and disinfecting mattresses. All cleaning and disinfecting solutions must be prepared according to the manufacturers' instructions. Follow standard institutional wipe down and infection control procedures. **DO NOT** open the zippered mattress top cover during cleaning, infection control procedures or during use. **DO NOT** launder or autoclave mattress covers.

Cover Material:

• Polyfiber/ Polyurethane Stretch

Wipe down the mattress unit with a damp cloth pre-soaked with warm water containing a mild detergent. Approved intermediate level disinfectants may be used according to the cover material. Please refer to the chart below for recommended disinfectants. After cleaning and disinfection, please avoid dust and proximity to dusty areas. All parts should be air dried thoroughly before use.

	Polyfiber/Polyurethane Stretch Cover	Vinyl Bottom Cover
Recommended	EPA-Approved Quaternary Ammonium Compound	Phenolic Disinfectants
Acceptable	Chlorine Bleach (1:10)	Quaternary, Chlorine Bleach (1:10)
Not Recommended	Phenolic Disinfectants	Isopropyl Alcohol

CAUTION: After cleaning, dry the mattress without direct exposure of sunlight.

1. Remove the bedding.
2. Clean the top and bottom covers with organic cleaner or soapy water.
3. Use a standard cleaning solution, described above, to clean the top and bottom covers.

3 Maintenance

3.1 Inspection

The internal and external parts of the mattress should be visually inspected every 3 months and in between patients for damage and wear. Remove the mattress from service and replace all damaged parts.

As an element of the maintenance, perform a compression test as per Section 3.2.

Inspection Record

INSPECTION DATE	INSPECTOR	COMPRESSION TEST RESULT	DAMAGE FOUND (Y/N)	REMOVED FROM SERVICE DATE	REPAIR DATE	RETURNED TO SERVICE DATE

3.2 Compression Test

ENTRAPMENT WARNING

Proper patient assessment and monitoring, as well as proper maintenance and use of mattresses is required to reduce the risk of entrapment. Visit the FDA website at <http://www.fda.gov> to learn about the risks of entrapment. Refer to the owners' manuals for beds and rails for additional product and safety information. Mattress MUST fit bed frame and bed rails snugly to reduce the risk of entrapment.

1. Remove the bedding from the mattress.
2. Measure the thickness of the mattress (mattress depth) in an uncompressed area.
3. Lay a straight edge over the compressed area of the mattress.
4. Measure the depth of the area of maximum compression (compression depth).
5. Use the following formula to determine the amount of compression (compression percentage):

$$\frac{\text{compression depth}}{\text{mattress depth}} \times 100 = \text{compression percentage}$$


FIGURE 1 Compression Test

6. Record the compression percentage calculated.

Note Compression percentages greater than 25% may be considered defective. Proof of this data point will be needed for any warranty claims.*

4 Limited Warranty

i PLEASE NOTE: THE WARRANTY BELOW HAS BEEN DRAFTED TO COMPLY WITH FEDERAL LAW APPLICABLE TO PRODUCTS MANUFACTURED AFTER JULY 4, 1975.

This warranty is extended only to the original purchaser/user of our products.

This warranty gives you specific legal rights and you may also have other legal rights which vary from state to state.

Drive warrants the mattress and cover when purchased new and unused to be free from defects in materials, workmanship, and assembly for the following time periods:

PRODUCT	ITEM #	FOAM	COVER
Gravity 6 Pressure Redistribution Mattress	15670	3 Years	3 Year

If within such warranty period any such product shall be proven to be defective, such product shall be repaired or replaced with new or refurbished product components, or a refund corresponding to the proper service time shall be offered, at Drive's option. This warranty does not include any labor or shipping charges incurred in replacement part installation or repair of any such product. Moreover, the warranty period of the repaired or replaced product shall terminate with the termination of the warranty period for the original product component.

For warranty service, please contact the dealer from whom you purchased your Drive product. In the event you do not receive satisfactory warranty service, please write directly to Drive at the address on the back cover. Provide the dealer's name, address, model number, and the date of purchase, indicate nature of the defect and, if the product has a lot number, indicate that as well.

Drive Corporation will issue a return authorization. The defective unit or parts shall be cleaned using disinfectant prior to return and must be returned for warranty inspection using the serial number, when applicable, as identification within thirty days of return authorization date. DO NOT return products to our factory without our prior consent. C.O.D. shipments will be refused; please prepay shipping charges. LIMITATIONS AND EXCLUSIONS: THE FOREGOING WARRANTY SHALL NOT APPLY TO LOT NUMBERED PRODUCTS IF THE LOT NUMBER HAS BEEN REMOVED OR DEFACED; CUSTOM SIZING OR ANY CUSTOMIZATION MADE TO THE FOAM AND/OR COVER; PRODUCTS SUBJECTED TO NEGLIGENCE, ACCIDENT, IMPROPER OPERATION, MAINTENANCE OR STORAGE; DAMAGE CAUSED BY ANY EVENTS OUTSIDE THE CONTROL OF DRIVE INCLUDING, BUT NOT LIMITED TO: DAMAGE FROM CARELESSNESS, ACCIDENTAL DAMAGE, DAMAGE FROM SHARP OBJECTS, TEARS, NEGLIGENCE (WHETHER WILFUL OR NOT), TAMPERING, DAMAGE FROM FIRE, DAMAGE OCCURRING IN TRANSIT AFTER DELIVERY BY DRIVE OR EXCESSIVE WEAR AND TEAR, INCLUDING BUT NOT LIMITED TO THE FAILURE TO FOLLOW CLEANING INSTRUCTIONS, USE OF INAPPROPRIATE CLEANING COMPOUNDS AND FAILURE TO ENSURE THAT ANY PRODUCT IS THOROUGHLY DRIED BEFORE TURNING OR FLIPPING IT OVER; DAMAGE RESULTING FROM CONTAMINATION FOLLOWING DAMAGE TO THE COVER AND FAILURE TO REPLACE IT, CONTAMINATION DUE TO INCORRECT STORAGE OR FAILURE TO INSPECT THE INTERIOR OF THE MATTRESS AS STATED HEREIN, OR DAMAGE DUE TO UNUSUAL CONTAMINATION OR INFECTION, SPILLAGE OR STAINING, OR USE OF ABRASIVE, CORROSIVE OR ANY SUCH SIMILAR DAMAGING SUBSTANCE; DAMAGE DUE TO UNDERTAKING ANY ACTION OR PROCEDURE USING THE PRODUCT OUTSIDE ITS NORMAL AND INTENDED USE; DAMAGE RESULTING FROM FAILURE TO REPLACE DAMAGED PARTS PROMPTLY OR DAMAGE DUE TO INTERIOR PARTS BEING USED WITHOUT THE COVER, BEING USED WITH A DAMAGED COVER, OR WITH THE COVER BEING INCORRECTLY FITTED; DAMAGE WHEN ANY PARTS OF THE MATTRESS ARE USED AS PART OF, OR IN CONJUNCTION WITH, ANY OTHER PRODUCT OR PIECE OF EQUIPMENT EXCEPT AS SPECIFICALLY PROVIDED FOR BY DRIVE; DAMAGE NOT PROVEN TO THE SATISFACTION OF DRIVE TO BE DUE TO DEFECTIVE MATERIALS, WORKMANSHIP OR ASSEMBLY; OR PRODUCTS MODIFIED WITHOUT DRIVE'S EXPRESS WRITTEN CONSENT INCLUDING, BUT NOT LIMITED TO MODIFICATION THROUGH THE USE OF UNAUTHORIZED PARTS OR ATTACHMENTS OR PRODUCTS DAMAGED BY REASON OF REPAIRS MADE TO ANY COMPONENT WITHOUT THE SPECIFIC CONSENT OF DRIVE . IT IS A CONDITION OF THIS WARRANTY THAT THE FOLLOWING IS COMPLIED WITH:THE FOREGOING EXPRESS WARRANTY IS EXCLUSIVE AND IN LIEU OF ANY OTHER EXPRESS WARRANTIES WHATSOEVER, WHETHER EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE AND THE SOLE REMEDY FOR VIOLATIONS OF ANY WARRANTY WHATSOEVER, SHALL BE LIMITED TO REPAIR, REPLACEMENT, OR REFUND OF THE DEFECTIVE PRODUCT PURSUANT TO THE TERMS CONTAINED HEREIN. THE APPLICATION OF ANY IMPLIED WARRANTY WHATSOEVER SHALL NOT EXTEND BEYOND THE DURATION OF THE EXPRESS WARRANTY PROVIDED HEREIN. DRIVE SHALL NOT BE LIABLE FOR ANY CONSEQUENTIAL OR INCIDENTAL DAMAGES WHATSOEVER AND DISCLAIMS ANY LIABILITY BEYOND THE INITIAL PURCHASE PRICE OF THE MATTRESS. THIS WARRANTY SHALL BE EXTENDED TO COMPLY WITH STATE/PROVINCIAL LAWS AND REQUIREMENTS.

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© 2022 Medical Depot, Inc. All rights reserved. Drive is a trademark of Medical Depot, Inc. Port Washington N.Y. 11050 USA Made in Mexico