# Harmony

**True Low Air Loss Tri-Therapy Mattress Replacement System**

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**Model Nos.: 14200, 14200-42 and 14200-48**

**PLEASE READ ALL INSTRUCTIONS BEFORE USE.**

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# IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE OPERATING THIS DEVICE.

## NOTE, CAUTION AND WARNING STATEMENTS DEFINED:

**NOTE** - Indicates helpful tips for using the system.

**CAUTION** - Indicates correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

**WARNING** - Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

**DANGER** - Indicates death or serious injury will result if proper precautions are not taken.

### **DANGER - To reduce the risk of electrocution:**

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

### **WARNING - To reduce the risk of burns, electrocution, fire, or injury to persons:**

1. Evaluate patients for entrapment risk according to facility protocol and monitor patients appropriately.
2. Use this product according to the manual and medical professionals' instructions.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accidents may result from a child swallowing a small part detached from the device.
4. Use this product only for its intended use as described in this manual.
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.
6. Keep the cord away from heated surfaces.
7. Never block the air openings of this product or place it on a soft surface, where their openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
8. Never drop or insert any object into any opening or hose.
9. Do not modify this equipment without authorization of the manufacturer.

10. Mattress covers have passed skin sensitization and skin irritation tests. However, If you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
11. Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.
12. Connect this product to a properly grounded outlet only. See Grounding Instructions.

**CAUTION -**

1. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (11 ft.) between devices or turn off the mobile phone.

**SAVE THESE INSTRUCTIONS**

**GROUNDING INSTRUCTIONS**

**DANGER** - Improper use of the grounding plug can result in a risk of electric shock.

If repair or replacement of the cord or plug is necessary, do not connect the grounding wire to either flat blade terminal. The wire with insulation having an outer surface that is green with or without yellow stripes is the grounding wire.

**NOTE** - If the repair or replacement of the cord is necessary, please contact a qualified electrician or serviceman. To reduce the risk of electric shock, do not modify the cord or plug in any way.

Check with a qualified electrician or serviceman if the grounding instructions are not completely understood, or if in doubt as to whether the product is properly grounded.

**SYMBOLS**

<b>IP20</b>	Protected against solid foreign objects of 12,5 mm and greater; No Protection against vertically falling water drops, Keep Dry! ( Only for 120V System)
<b>IP21</b> 	Protected against solid foreign objects of 12,5 mm and greater; Protection against vertically falling water drops ( Only for 230V System)
	Temperature limitation/temperature range
	Consult operating instructions for use.
	"BF" symbol, indicate this product is according to the degree of protecting against electric shock for type BF equipment
	Grounding terminal
	Attention, read the accompanying information carefully.
	Do Not Bleach
	Do Not Iron
	Tumble Dry, Normal, Low Heat
	Dry clean, Any Solvent Except Trichloroethylene
	Machine wash, regular / normal, 95 degrees C (203 degrees F)
	Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an applicable collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

# 1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

## 1.1 GENERAL INFORMATION

The system is a Tri-Therapy mattress replacement system including alternating pressure, pulsation and true low air loss. These therapies are intended for the prevention and treatment of all stages of pressure injuries.

The system has been tested and successfully approved with the following standards:

- EN 60601-1
- EN 60601-1-2
- EN 55011 Class B
- IEC61000-3-2
- IEC 61000-3-3

**FOR US AND CANADA ONLY**

ETL CLASSIFIED

**Intertek**

3176994

CONFORMS TO AAMI  
STD. ES60601-1  
CERTIFIED TO CSA  
STD. C22.2 NO. 60601-1

The product was tested with medical equipment and meets AAMI STD. ES60601-1 & CSASTD. C22.2 NO. 60601-1. preventing electric shock, fire and risk of physical injury.

### EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation distances between the equipments.
- Connect the equipment into an outlet on a circuit different from that to which the other device is connected.
- Consult the manufacturer or field service technician for help.

## 1.2 Intended Use

This product is intended:

- to help reduce the incidence of pressure injuries while optimizing patient comfort
- for use primarily in acute or long term health care facilities
- for pain management


NOTE: Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

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## 2. Product Description


### 2.1 Pump and Mattress System

1. Quick Connector for CPR
2. Control Unit
3. Mattress
4. Mattress Foam Pocket
5. Cable Management Loops


## 2.2 Control Unit


Front

1. Quick Connector for CPR
2. Front Panel


Rear

1. Power Switch
2. Power Socket
3. Fuses
4. Integrated Handle
5. Mounting Bracket
6. Pump Bumper
7. Air Filter

## 2.2 Control Panel

## 1. Comfort Weight Settings Buttons

The “Comfort Setting” controls the air pressure output based on the patient’s weight. When patient’s weight setting is increased (  ), the output pressure will increase. When the patient’s weight setting is decreased (  ), the output of air pressure will decrease. When a patient’s condition has significantly changed, reassess the comfort weight setting level. When a patient’s condition has significantly changed, reassess appropriateness of product and comfort setting level.

**NOTE: Maximum weight: 14200= 440 lbs**

**14200-42= 700 lbs**

**14200-48= 990 lbs**

## 2. Therapy Modes Button

### A. Static

Redistribute body mass over a greater surface area at a constant air pressure based on your Comfort Weight setting.

### B. Alternating Pressure

1-in-2 alternating cell cycle achieves periodic pressure relief. There are four selectable cycle times available: 5, 10, 15, 20 minutes. Caregivers can select one of the four cycle times based upon patient comfort and desired outcome.

### C. Pulsation

Encourage lymph and blood flow for increased oxygenation. The entire mattress pulsates between a decreased and increased pressure level in all cells every 15 seconds.

### D. Seat Inflate

The Seat Inflation features additional support to the patient during head raised position (Fowler), ingress, and egress for the prevention of bottoming out.

**NOTE:** When head of bed is elevated 30 degrees or higher, it is recommended the caregiver engage Seat Inflate to prevent possible bottoming out.

**NOTE:** Every time the mattress is initialized (inflated), it will automatically go to “Auto Firm” mode to hasten inflation. Once the system is ready to use, the system will automatically begin in static mode or revert to the previously set therapy mode.

## 3. On/off

Press to turn the unit On or Off. Press Off to turn Off/Standby the unit.

**NOTE:** The power switch on the side of pump must be turned on.

#### 4. Panel Lock

Should the panel remain untouched for 5 minutes, the panel lock feature will Lock the panel and a green LED light will display. This feature helps prevent accidental or unintentional setting changes. To unlock, press Panel Lock button for 3 seconds.

#### 5. Cycle Time

During alternating pressure therapy mode there are four selectable cycle times: 5, 10, 15, 20 minutes. By pressing the Cycle Time button, User can select one of the four cycle times based upon patient comfort and desired outcome.

#### 6. Alarm Mute

Press alarm mute button to temporary suspend the alarms for low pressure, power failure or system failure. Should the problem continue, the alarms will resound within three minutes to notify the caregiver until resolved.

##### A. Low Pressure Indicator

Should the Quick connector become disconnected, the lower pressure LED indicator will light up and the low pressure alarm will activate. This alarm may sound within 20 seconds. Once the low pressure problem is resolved the Control Unit will resume operation in the previously set mode. NOTE: The low pressure indicator will not be detected if air cells of mattress are purposely disconnected.

##### B. Power Failure Alarm

Should there be a power failure outage, the Power Failure LED indicator will light up and the Power Failure alarm will activate immediately. Press the alarm mute button to deactivate.

##### C. System Failure

Should there be a System Failure issue, the System Failure LED indicator will light up. This may be caused by a motor related issue. Please consult the manufacturer or field service technician for help.


#### 7. AutoFirm/ Nursing Key

The system will go into AutoFirm mode automatically when the power switch is turned on. This insures the Control Unit reaches its maximum operating pressure. Once the maximum pressure level is reached, the Control Unit will automatically switch to the previously selected comfort level in **STATIC** mode, or the user can press the **THERAPY** modes button to select another therapy. User can also use this function as full mattress inflation during patient ingress/egress or normal nursing procedure for better support.

### 3. Installation

Unpack the box to inspect for any damage, which may have occurred during shipment. If there are any damages, please contact your provider immediately.


#### 3.1 Control Unit & Mattress Installation

1. Place the mattress on top of the bed frame. Please note for the foot end. There are securing straps on the base of the mattress. Secure the mattress firmly by fixing the straps to the bed frame, ensure that moving sections of the bed are still free to move.


**WARNING:** The Harmony mattress is a mattress replacement system and must not be installed on top of an existing mattress.

2. Hang the control unit onto bed rail (foot-end), and adjust hangers to best upright position of the Control Unit.

**NOTE:** Do not place the Control Unit on the floor.

3. Connect the Quick Connector from air mattress to the Control Unit. Make sure the connector is in the right position as per the diagram on the left. When a “click” sound is felt or heard, the connection is completed and secured.

**NOTE:** The Control Unit is only operable when the Quick Connector is connected to the system.

**NOTE:** Check and ensure the air hoses are not kinked or tucked under mattress.


4. Plug the power cord into electrical outlet.

**NOTE:**

Make sure the Control Unit is suitable for the local power voltage.


The plug is also served to disconnect the device. Do not position the equipment so that it is difficult to disconnect the device.

**CAUTION:** The Control Unit can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose.

5. Then turn the main power switch to ON position.

6. After installation, the extra length of the power cord, if any, should be neatly arranged through the cable management loops on the left side of the mattress to avoid a tripping hazard.

**NOTE:** Do not connect any cable to the reserved signal port next to the quick connector. The cap should always remain on.


**Grounding:**

See Grounding Instructions on Page 2.


## 4. Operation

NOTE: Always read the operating instructions before use.

### 4.1 GENERAL OPERATION


1. Switch on the main power switch found on the left side of the Control Unit while facing the User. Press (⏻) on the control panel to turn on the power.
2. The system will automatically go into “Auto Firm” mode for a few minutes of inflation.
3. When the initial inflation (Auto Firm process) is completed, the system will automatically enter into Static mode
4. Using the Comfort Weight Setting buttons, adjust according to the weight and height of the patient, adjust the pressure setting to the most suitable level without bottoming out.
5. Using the Therapy Mode button choose the desired therapeutic mode.


**NOTE:** Seat inflate mode should be used for ingress and egress at all times to prevent bottoming out.

### 4.2 CPR

When CPR needs to be performed, click on the release button on Quick Connector and quickly detach the connector from the system to release air.

## 5. Cleaning

It is important to follow the cleaning procedures before using the system.

Wipe down the Control Unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the Control Unit. If your doctor or medical facilities have other special cleaning instruction, please follow the professional instruction.

**CAUTION:** Do not immerse or soak Control Unit.

**WARNING:** Do not remove the housing of the Control Unit to avoid electrical shock. All disassembly or repair should be done by professional technicians.

**CAUTION:** The Control Unit does not need oil lubrication, please do not disassemble the system.

<p><b>Cover Material:</b></p> <ul style="list-style-type: none"> <li>• Silver<sup>®</sup> Ion</li> <li>• Polyfiber/ Polyurethane Stretch</li> </ul>	
---	--

Wipe down the mattress unit with a damp cloth pre-soaked with warm water containing a mild detergent. Approved intermediate level disinfectants may be used according to the cover material. Please refer to the chart below for recommended disinfectants.

The mattress top cover can be completely removed for laundry with water temperature up to 95°C; however, it is recommended that the user still check with local policy to determine the time/temperature ratio required to achieve thermal disinfection. After cleaning, please avoid dust and proximity to dusty areas. All parts should be air dried thoroughly before use.

	Polyurethane Stretch Top Cover	Vinyl Bottom Cover
Recommended	Quaternary	Phenolics
Acceptable	Chlorine Bleach (1:10)	Quaternary, Chlorine Bleach (1:10)
Not Recommended	Phenolics	Quat/Isopropyl

**CAUTION:** After cleaning, dry the mattress without direct exposure of sunlight.

## 6. Storage

1. To quickly deflate the mattress for storage, click on the release button on the Quick Connector
2. Roll from the foot end towards the head end with CPR valve open, and make sure the tubing is not kinked.
3. Foot-end strap can then be stretched around the rolled mattress to prevent unrolling. Fasten the buckle strap to secure the packed mattress.
4. The power cord can be wrapped around the Control Unit bumper on the back pack the Control Unit with protective package.
5. Place the whole system into the carrying bag.

**NOTE:** Do not kink, crease or stack the mattresses and do not store the system in direct sunlight, high temperature or moist area.


## 7. Maintenance

### 7.1 General

1. Check main power cord and plug for abrasions or excessive wear.
2. Check mattress cover for signs of wear or damage. Ensure mattress cover is securely on and tubes are connected together correctly.
3. Check the air hoses for any kinks or breaks. For replacement, please contact your local provider.


### 7.2 Fuse replacement

#### External Fuse replacement:

1. Disconnect the plug from the main power source when a blown fuse is suspected.
2. Use a proper tool to open the fuse holder.
3. Insert a new fuse of the correct rating in, and replace the fuse holder back to correct position. The fuse should be rated as 5A / 250V (for 120V system) and T3.15A /250V (for 230V system) type.

### 7.3 Air Filter Replacement

1. Replace the air filter located at the side of the Control Unit.
2. The filter is reusable and can be washed gently with a mild detergent and water. Air dry the filter before use.
3. Check and replace air filter regularly if environment is dirty

### 7.4 Rechargeable Battery

1. The rechargeable battery is designed to support power failure alarm. To check if the rechargeable battery has been drained out, unplug the power cord and see if the Power Failure indicator will light up along with buzzer and last for a few minutes.
2. If the Power Failure Alarm is unable to work or the battery might need to be replaced (approximate life expectancy 6 months or 300 cycles), please contact your dealer or notify the technician for replacement.

## 8. Troubleshooting

Problem	Solution
Power is not ON	Check if the plug is connected to control unit Check for a blown fuse.
Low Pressure Alarm is on	Check if the Quick Connector is tightly secured. Check if all tubing connections along mattress are secured. Check if the air hoses are kinked or broken
Power Failure Alarm is on	Check if the power is suddenly shut down. Check if the power cord is connected properly.
Patient is bottoming out	Pressure setting might be inadequate for the patient, adjust comfort range 1 to 2 levels higher and wait for a few minutes for best comfort.
Mattress form is loose	Check if all the snap buttons or straps of mattress are all securely fastened. Check if the mattress is fixed to the bed frame by straps.
No air produced from some air outlets of the air tube connector	This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.

If the above information does not solve your problems, please contact your local agent directly. They might require a technician to take care the problem.

## 9. EXPECTED SERVICE LIFE

This product is intended to offer safe and reliable therapy when installed and operated according to the instructions provided by Drive Medical. Drive Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device functionality and indications on products. Otherwise, service and inspection of the devices generally should not be required.

# 10. Technical Description

## SPECIFICATIONS:

Pump				
Power Supply (Note: See rating label on the product)	AC 110-120V 60Hz, 4.5A (for 120V system)			
Fuse Rating	T5AL 250V(for 120V system)			
Dimension (L x W x H)	34 x 16.3 x 23.5 (cm) or 13.4" x 6.4" x 9.3" (in)			
Weight	5 kg or 11.02 lb			
Classification	Class I, Type BF, IP21 (for 230V System ), IP20 (for 120V System) Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)			
Environmental Information				
Environment	Temperature	Operation: 10° C to 33° C (50° F to 91° F) Storage: -15° C to 50° C (5° F to 122° F) Shipping: -15° C to 70° C (5° F to 158° F)		
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping:10 % to 90% non-condensing		
	Operating Altitude	2000 m (Maximum)		
Mattress				
Model	8" Mattress 2" foam pocket is provided			
Size of Mattress	36" mattress	42" mattress	48" mattress	
Dimension (L x W x H)	78.7" x 35.4" x 8" 200 x 90 x 20.3 cm	78.7" x 42" x 8" 200 x 107 x 20.3 cm	78.7" x 48" x 8" 200 x 122 x 20.3 cm	
Weight	16.5 lb / 7.5 kg	18.8 lb / 8.5 kg	21.0 lb / 9.5 kg	
Maximum Weight Capacity	Static mode	440 lb / 200 kg	700 lb / 320 kg	990 lb / 450 kg
	Therapy mode	400 lb / 180 kg	650 lb / 300 kg	920 lb / 420 kg

- NOTE:**
1. Consult Drive Medical for further technical documents.
  2. The specification is also suitable for other areas operating with same power supply.
  3. Mattress dimension and weight is measured without foam cushion
  4. The manufacturer reserves the right to modify the specification without notice.

# Appendix A: EMC Information

## Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	
<b>WARNING</b> 1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used. 2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. 3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PUMP, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.		

## Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC Standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional Healthcare Facility Environment	Home Healthcare Environment		
Electrostatic Discharge (ESD) IEC61000-4-2	±8kV contact ±5kV air		±8kV contact ±5kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/output line		±2kV for power supply line ±1kV for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s)	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of atypical commercial or hospital environment.

Basic EMC Standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional Healthcare Facility Environment	Home Healthcare Environment		
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25/30 cycles <5 % UT (>95 % dip in UT) for 250/300 cycles		120V	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz	6Vrms	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM Fields IEC61000-4-3	3 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz	3V/m	<p>Recommended separation distance  <math>d=2.0 \cdot P^{0.5}</math> 150kHz to 80MHz  <math>d=2.0 \cdot P^{0.5}</math> 80MHz to 800MHz  <math>d=4.0 \cdot P^{0.5}</math> 800 MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency ranged.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
<p>NOTE 1: <math>U_r</math> is the a.c. mains voltage prior to the application of the test level</p> <p>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p>				

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and this device:**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 2.0\sqrt{P}$	80 MHz to 800 MHz $d = 2.0\sqrt{P}$	800 MHz to 2,7 GHz $d = 4.0\sqrt{P}$
0.01	0.2	0.2	0.4
0.1	0.63	0.63	1.26
1	2	1	4
10	6.3	6.3	12.6
100	20	20	40

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

**Warranty**

**Drive Medical**  
**99 Seaview Blvd.**  
**Port Washington, NY 11050**

**Phone: 1-877-224-0946**  
**www.drivemedical.com**

Your Drive branded product is warrantied to be free of defects in materials and workmanship for 2 years from the date of purchase for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Limited 2 Year Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.

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