DeVilbiss® PulmoNeb® LT Compressor Nebulizer System Instruction Guide

Model 3655LT

CAUTION– Federal (U.S.A.) law restricts this device to sale by, or on the order of a physician.

Not made with natural rubber latex.

Conforms to AAMI Std. ES60601-1, IEC Std. 60601-1-6

Made in China

IEC SYMBOLS

	It is mandatory to read Operating Instructions prior to use.		Type BF Applied Part
	Attention Consult Instruction Guide		Double Insulated
	Electric Shock Hazard		Catalog/Model Number
	Alternating Current		Serial Number
	Off On		Lot Number
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU - Waste Electrical and Electronic Equipment (WEEE)		

READ ALL INSTRUCTIONS BEFORE USING.
Accompanying documents specify special skills, training, and knowledge required of operator or responsible organization and environmental restrictions on locations of use.

DANGER

To reduce the risk of electrocution:

- Always unplug this product immediately after using.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink or other water source.
- Do not reach for a product that has fallen into water. Unplug immediately.
- DeVilbiss compressor/nebulizer must be operated on the specified power source to avoid risk of electric shock and damage to the compressor.
- Do not submerge in water; doing so will result in compressor damage.
- Electric shock hazard. Do not remove compressor cabinet. All disassembly and maintenance must be done by a qualified DeVilbiss provider.

WARNING

To reduce the risk of burns, electrocution, fire or injury to persons:

- The product should never be left unattended when plugged in.

- Close supervision is necessary when this product is used by, on, or near children or physically incapacitated individuals.
- Use this product only for its intended use as described in this guide. Do not use attachments not recommended by the manufacturer.
- Never operate this product if:
 - It has a damaged power cord or plug.
 - It is not working properly.
 - It has been dropped or damaged.
 - It has been dropped into water.Return the product to an authorized DeVilbiss service center for examination and repair.
- Keep the power cord away from heated surfaces.
- Never block the air openings of the product or place it on a soft surface, such as a bed or couch, where the air openings may be blocked. Keep the air openings free of debris.**
- Never use while drowsy or asleep.
- Never drop or insert any object into any opening or hose.
- Do not use outdoors. This product is for household use only.
- Do not use in an oxygen enriched environment.
- Do not use this product with a DC to AC converter, or with any AC voltage and/or frequency other than specified.
- Do not operate unit in a backpack or carrying bag.
- Use caution when plugging in and unplugging device. Do not allow jewelry or other metal items to contact plug.
- This product contains small parts that may represent a choking hazard.**

- To prevent possible risk of infection from contaminated medication, cleaning of the nebulizer is recommended after each aerosol treatment. Daily disinfecting is recommended.
- To prevent possible risk of infection from contaminated cleaning solutions, always prepare fresh solution for each cleaning cycle and discard solution after each use.


NOTE–This unit employs a temperature and/or current limiting device that will shut off the entire product should a malfunction occur or in the event of product misuse, abuse, or abnormal operation.

NOTE–The exterior surfaces of the unit and the air outlet connector may become slightly warm to the touch during use. Warm air can also be felt exiting the vents on the sides of the case of the unit. This is normal.

SAVE THESE INSTRUCTIONS STATEMENT OF INTENDED USE

The DeVilbiss compressor nebulizer includes an AC powered air compressor that provides a source of compressed air for home health care use. The compressor is used in conjunction with a jet (pneumatic) nebulizer to convert certain inhalable drugs into an aerosol form for inhalation by a patient. The device can be used with adult or pediatric patients.

IMPORTANT PARTS OF YOUR DEVILBISS 3655LT COMPRESSOR AND NEBULIZER

COMPRESSOR

- Power switch
- Carrying handle
- Nebulizer holder
- Filter (non-washable)
- Air-outlet connector

NEBULIZER *

- Mouthpiece
- T-piece
- Corrugated Tubing
- Cap
- Baffle
- Medication Cup
- Nebulizer Air-Inlet Connector
- Tubing

* **NOTE**–Depending on model purchased, your unit may include a disposable or reusable nebulizer, or both.

HOW TO OPERATE YOUR COMPRESSOR

NOTE–Before initial operation, clean nebulizer.

- Ensure the power switch is in the “O” off position.
- Place compressor on a level, sturdy surface. Never block air vents. Never place compressor on a soft surface, such as a patient’s lap, or a bed or couch, where the air openings may be blocked.
- Plug power cord into wall outlet.
- Wash hands.
- Connect one end of the nebulizer tubing to the compressor air-outlet connector.

NOTE–During high-humidity weather conditions, condensation (water build up) can occur in the nebulizer tubing.

- Assemble nebulizer as shown. If using an aerosol mask, insert the bottom part of the mask directly into the top of the nebulizer cap.
- Press the power switch to start compressor.
- Place the mouthpiece between teeth. With mouth closed, inhale deeply and slowly through mouth as aerosol begins to flow, then exhale slowly through the mouthpiece. If treatment needs to be interrupted, simply press power switch and place the nebulizer into nebulizer holder.
- If an aerosol mask is used, place the mask over mouth and nose. As aerosol begins to flow, inhale deeply and slowly through mouth, then exhale slowly.
- When treatment is complete, turn unit off. Unplug unit from wall outlet.

NEBULIZER CLEANING

- Disassemble nebulizer.
- Wash all items, except tubing, in a hot water/dishwashing detergent solution. Rinse under hot tap water for 30 seconds to remove detergent residue. Allow to air dry.

Disinfect Daily:

- Using a clean bowl, soak items in three parts hot water (131°F-149°F/55°C- 65°C) to one part white vinegar (>=5% acetic acid concentration) for 30 minutes OR use a medical bacterial-germicide disinfectant. Be sure to follow manufacturer’s instructions carefully.
- With clean hands, remove items, rinse under hot tap water, and air dry. Store in a zip-lock bag.
- Keep the outer surface of the tubing dust-free by wiping regularly. Nebulizer tubing does not have to be washed because only filtered air passes through it.

NOTE–Replace your DeVilbiss primary reusable nebulizer every 6 months. Disposable nebulizer (often used as a backup) should be replaced at least every 15 days. To ensure consistent operating specifications, replace with DeVilbiss nebulizer.

COMPRESSOR CLEANING

- With power switch in the “O” off position, unplug power cord from wall outlet.
- Wipe outside of the compressor cabinet with a clean damp cloth every few days to keep dust-free.

Filter Change

- Filter should be changed every 6 months or sooner if filter discolors.
- Remove filter by grasping it firmly and pulling out of the unit. Discard old filter.
- Replace with a new DeVilbiss filter (part #3655LT-601). Additional filters should be purchased from your DeVilbiss provider.

CAUTION–Do not reuse or wash filter or substitute any other material for a DeVilbiss air-inlet filter; it will result in compressor damage.

MAINTENANCE

All maintenance must be performed by a qualified DeVilbiss provider.

This unit is not Field serviceable; it will be replaced, not repaired. Please refer to warranty statement.

SPECIFICATIONS

Size	5.75” W x 4.0” H x 6.0” D (14.6 cm x 10.2 cm x 15.2 cm)
Weight	2.75 lbs. (1.25 kg)
Electrical Requirements	1.2 amps, 120 V~ 60 Hz
Sound Level	56 dBA
Power Consumption	110 VA maximum
Maximum Compressor Pressure	30 psig or greater (206 kPa)*
Compressor Free Air Flow	12.7*
Operating Temperature Range	+40° to +104°F (+5 to +40°C)
Operating Humidity	up to 95% non-condensing
Storage/Transport Temperature Range	–40° to +158°F (–40° to +70°C)
Storage/Transport Humidity	up to 95% non-condensing
Supplied Nebulizer Capacity	6 ml (cc)
Supplied Nebulizer Nebulization Rate	0.15 ml/min or greater*
Supplied Nebulizer MMAD	5 microns or less**

*Conditions may vary based on altitude above sea level, barometric pressure, and temperature.

**Value determined with minimum flow rate through the nebulizer (20°C, sea level).

There are no hazardous materials used in the construction of this unit, therefore, there are no known risks associated with disposal of the equipment.

This product complies with the following electromagnetic compatibility standard: IEC 60601-1-2: 2014 3rd Edition.

Equipment Classifications:

Protection against electric shock: Class II
Degree of protection against electric shock:.....Type BF
Equipment not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.
Mode of operation:.....Continuous operation

FIVE-YEAR LIMITED WARRANTY ON COMPRESSOR PARTS

All DeVilbiss Compressor models are warranted to be free from defective workmanship and material for a period of five years from the date of purchase. The unit will be replaced if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your authorized DeVilbiss provider or DeVilbiss Service Department at 800-338-1988 (814-443-4881) for instructions. **NOTE**–This warranty does not cover providing a loaner compressor, compensating for costs incurred in compressor rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

DEVILBISS GUIDANCE AND MANUFACTURER’S DECLARATION

WARNING

Medical Electrical Equipment needs special precautions regarding Electromagnetic Compatibility (EMC) and needs to be installed and put into service according to the EMC information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

NOTE–The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

Guidance and Manufacturer’s Declaration – Emissions All Equipment and Systems			
This device is intended for use in the electromagnetic environment specified below. The customer or user of this device should assure that it is used in such an environment.			
Emissions Test	Compliance	Electromagnetic Enforcement – Guidance	
RF Emissions CISPR 11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class B Radiated and Conducted Emissions	This device is suitable for use in all establishments including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonics IEC 61000-3-2	Class A		
Flicker IEC 61000-3-3	Complies		
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8kV contact ±16kV air	±8kV contact ±16kV air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
Electrical Fast Transient/burst IEC 61000-4-4	±2kV on AC Mains	±2kV on AC Mains	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95% Dip for 0.5 Cycle 30% Dip for 25 Cycles >95% Dip for 5 Seconds	>95% Dip for 0.5 Cycle 30% Dip for 25 Cycles >95% Dip for 5 Seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply.
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be that of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms from 150 kHz to 80 MHz	V1 = 3 Vrms	Portable and mobile RF communications equipment should be separated from the device by no less than the recommended separation distances calculated/listed below: D=(1.2)√P
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	E1 = 3V/m	D=(1.2)√P 80 to 800 MHz D=(2.3)√P 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and this device. This device and system are NOT Life-Supporting

This device is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Maximum Output Power (Watts)	Recommended Separation Distances for the device (meters)		
	150 kHz to 80 MHz D=(1.2)√P	80 to 800MHz D=(1.2)√P	800 MHz to 2.5 GHz D=(2.3)√P
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.30
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Manufactured For/ Fabricado para:

DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501-2125 • USA
800-338-1988 • 814-443-4881 • www.devilbisshealthcare.com

DeVilbiss® and PulmoNeb® are registered trademarks of DeVilbiss Healthcare.

© 2016 DeVilbiss Healthcare LLC. 10.16 All Rights Reserved.

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