Compressor Nebulizer System

Model MQ5800, MQ5900, MQ5900P

Instruction Guide

Sistema de Compresor
Nebulizador
Modelo MQ5800, MQ5900, MQ5900P

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Guide d'Instructions

Caution– Federal (U.S.A.) law restricts this device to sale by, or on the order of a physician.

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Introduction.....6

 Statement Of Intended Use.....6

 What Are Compressor Nebulizer Systems Used For?.....6

 How Does A Compressor Nebulizer Work?.....6

 How Do I Take A Treatment?.....7

**Product Identification
and Available Accessories.....8**

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OFF

IP 21

IP21 — Protected against foreign objects equal to greater than 12.5 mm in diameter and against vertically falling drops of water

Alternating Current

Type BF Applied Part

Rx only: Federal (U.S.) Law restricts this device to sale by or on the order of a physician

WEEE
Waste Electrical
and Electronic
Equipment

Class II Equipment

Attention, consult
instruction guide.

Refer to instruction
manual / booklet

Important Safeguards

When using electrical products, especially when children are present, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

ANGER

Urgent safety information for hazards that will cause serious injury or death.

ARNING

Important safety information for hazards that might cause serious injury.

UTION

Information for preventing damage to the product.

OTE

Information to which you should pay special attention.

Read All Instructions Before Using

To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

**DANGER!**

ARNING:

If any of the motor housing parts should become detached or broken exposing the motor or any other electrical component operation should be discontinued immediately to avoid personal injury or further product damage.

To reduce the risk of burns, electrocution, fire or injury to persons:

**WARNING!**

1. The product should never be left unattended when plugged in.
2. Close supervision is necessary when this product is used by, on, or near children or disabled individuals
3. Use this product only for its intended use as described in this guide. Do not use attachments not recommended by the manufacturer.
4. Never operate this product if:
 - a. It has a damaged power cord or plug.
 - b. It is not working properly.
 - c. It has been dropped or damaged.
 - d. It has been dropped into water.
5. Keep the power cord away from heated surfaces.
6. Never allow the power cord to be pulled, jerked, strained, twisted or severely bent

especially at the plug connections. Damage will occur at the high flex point of entry to the device, causing it to rupture and short.

7. Never block the air openings of the product or place it on a soft surface, such as bed or couch, where the air openings may be blocked. Keep the air openings free of lint, hair, and the like.
8. Never use while drowsy or asleep.
9. Never drop or insert any object into any opening or hose.
10. Do not use outdoors. This product is for household use only.
11. DO NOT use this product in or near explosive environments, or where aerosol (spray) products are being used.
12. DO NOT pump combustible liquids or vapors with the product or use in or near an area where flammable or explosive liquids or vapors may exist.
13. DO NOT operate this product in an oxygen enriched environment, i.e. oxygen tent or oxygen hood.
14. Use only in well ventilated areas.
15. Do not use this product with a DC to AC converter, or with any AC voltage and/or frequency other than specified.
16. Do not operate unit in a backpack or carrying bag.
17. Care should be taken when this product is used in the vicinity of other electronic equipment concerning the possible interference between them. They may need to be moved apart.
18. Use caution when plugging in and unplugging device. Do not allow jewelry or other metal items to contact plug.
19. Warning: This product contains small parts that may represent a choking hazard.
20. The device does not need to be calibrated.
21. Do not position the equipment so that it is difficult to operate.
22. To avoid strangulation and entanglement, keep power cord and air tubes out of reach of young children.
23. The air filter shall not be serviced or maintained while the compressor is in use.
24. Warning: No modification of this equipment is allowed.

NOTE This unit employs a temperature and/or current limiting device that will shut off the entire product should a malfunction occur or in the event of product misuse, abuse, or abnormal operation.

NOTE This is based on the thermal interrupt safety device that is integrated into the motor. As the unit runs it generates heat and gets warm. If the unit is run for over 30 minutes continuously the thermal interrupt may trip which would stop the user's treatment. Some models of this appliance have a polarized plug (one blade is wider than the other). As a safety feature, this plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician. Do not attempt to defeat this safety feature. If repair or replacement of the cord or plug is necessary, contact your qualified Drive provider.

NOTE The exterior surfaces of the unit and the air outlet connector may become slightly warm to the touch during use. Warm air can also be felt exiting the vents on the sides of the case of the unit. This is normal. Your doctor has prescribed a liquid medication to treat your respiratory condition. To best utilize this liquid medication, he has prescribed a Drive brand compressor nebulizer. Your Drive compressor nebulizer works to convert the medication to a high-quality mist of fine particles that penetrates deep into the lungs.

Be sure that you read and understand the information in this instruction guide. By following these simple instructions and the advice of your physician, your compressor will become an effective addition to your therapeutic routine.

Statement of Intended Use

The aerosol generator is for the production of compressed air to operate a nebulizer for the production of medical aerosol for respiratory disorders. The aerosol generator is used in conjunction with a pneumatic nebulizer which nebulizes specific drugs for inhalation by a patient and requires the order of a physician for medical use.

What are Compressor Nebulizer Systems Used For?

People with breathing problems, such as Asthma, Chronic Obstructive Pulmonary Disease (COPD) or Cystic Fibrosis are often prescribed medicines that minimize or control symptoms of their disease. Compressor Nebulizer Systems are a means of delivering these prescribed medications to the lungs. These medications and the compressor nebulizer systems should be used as directed by a physician to maximize benefits. Compressor Nebulizer systems can be prescribed for use everyday or periodically as symptoms present themselves, such as the springtime when Hay Fever or Allergies are normally present.

How Does a Compressor Nebulizer Work?

Most medicines for use with a compressor nebulizer system are created in a liquid form. This liquid medication is placed in the “cup” of the nebulizer (please refer to this Instruction Guide for more information on the proper use of a nebulizer). The compressor, which is connected to the nebulizer by flexible tubing, generates a specific pressure and flow of air that essentially “drives” the nebulizer. The pressurized air produced by the compressor initiates the nebulizer’s ability to break down the liquid medication into a very fine mist. This mist, often referred to as an aerosol cloud, is then propelled out the nebulizer and inhaled into the lungs where the medicine is readily absorbed. Please consult your prescription for details on the amount of liquid medication that should be placed in the nebulizer. If you should have any questions about your medication or its use, please contact your physician or pharmacist.

How Do I Take a Treatment?

A nebulizer is the optimal delivery method of medication to the lungs. However, some users may find it difficult to use, especially very young children. In this case, an aerosol mask may be used. To ensure that the proper mask is utilized, please consult your home health provider or physician.

The most important aspect of taking a treatment is the patient's breathing pattern while inhaling the medication mist.

It should be slow and deep, all attempts should be made to sit up straight and breathe in deep.

A general rule of thumb is to attempt to breathe in so deeply that it feels as though the lungs are completely full of air and inhaling any additional air is impossible.

If possible, hold your breath for 1-2 seconds unless directed otherwise by your physician and then slowly exhale through the nebulizer mouthpiece or mask. Holding your breath for 1-2 seconds after inhaling the medication allows for optimal absorption of medication in the lungs.

Repeat this breathing pattern for the entire treatment.

The length of a treatment can last anywhere from 5 to 20 minutes, depending on the amount of medication and the performance of your compressor and nebulizer system. When you cannot see any medication in the nebulizer cup and the nebulizer has stopped producing a mist, the treatment has been completed. Turn off your compressor and follow the cleaning instructions found in this instruction guide.

If you have any additional questions regarding your equipment, medication, or treatment, please consult your physician.

Product Identification

MODELS	DESCRIPTIONS
MQ5800	Compact Compressor Nebulizer with JetNeb Kit
MQ5900	Compact Compressor Nebulizer with JetNeb Kit and Neb Kit 500
MQ5900P	Compact Compressor Nebulizer with JetNeb Kit, Neb Kit 500 and JetNeb Pediatric Mask

Available Accessories

18260	18260 JetNeb Kit with mouthpiece, tubing (12/cs)
NEB KIT 500	Disposable Nebulizer Kit with mouthpiece, 7' tubing (50/cs)
MASK 001P	Pediatric Aerosol Mask (50/cs)
MQ92692	Replacement Filter (5/bag)

How to Operate Your Compressor

NOTE

Before initial operation, the nebulizer cup assembly should be cleaned following instructions described in the “Cleaning” section.

WARNING

Before connecting the power cord, make sure the I/O (ON/OFF) switch is in the O (OFF) position.

The plug is also served to disconnect the device.

Daily Use Operation

CAUTION

The Compressor is designed for intermittent use only. Do not operate it continuously for more than 30 minutes for a single use without turning it off and following a cooling period for least 30 minutes.

1. Before each use inspect the Compressor and nebulizer cup assembly for damage or wear, replace as needed.
2. Place the Compressor on a table or other flat stable surface. Be sure you can easily reach the controls when seated. Do not use this device on the floor.
3. With the power switch in the O (OFF) position, plug the power cord into an appropriate electrical wall outlet.
4. Connect one end of the tubing to the compressor air-outlet connector.
5. Assemble the nebulizer cup and add indicated medication to the nebulizer's cup before use.
6. Attach the other end of the tubing into the air-inlet connector found at the bottom of the nebulizer cup.
7. Turn on the Compressor by pressing the power switch to the I (ON) position and begin treatment.
8. If treatment needs to be interrupted, simply press power switch to O (OFF) position.
9. When the treatment is complete, turn off the compressor by pressing the power switch to O (OFF) position and unplug the unit from the electrical outlet.
10. Always unplug this device immediately after use.

Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. This device has no AP or APG protection not intended for use with oxygen or oxygen enriched environments.

11. DO NOT use this product in or near explosive environments, or where aerosol (spray) products are being used.
12. DO NOT pump combustible liquids or vapors with the product or use in or near an area where flammable or explosive liquids or vapors may exist.
13. DO NOT operate this product in an oxygen enriched environment, i.e. oxygen tent or oxygen hood.
14. Use only in well ventilated areas.

Safety Overload

NOTE Do not exceed 30 minutes of continuous operation.

NOTE This is based on the thermal interrupt safety device that is integrated into the motor. As the unit runs it generates heat and gets warm. The thermal interrupt will cut power to the unit if the temperature inside gets to an unsafe level. Once the unit cools off it will work normally. If the unit is run for over 30 minutes continuously the thermal interrupt may trip which would stop the user's treatment

1. The motor of this device has a built-in thermal overload protector. Should the motor overheat, the protector will automatically shut off the motor. Should this occur, turn the I/O (ON/OFF) switch to the O (OFF) position and allow the motor to cool down for approximately 30 minutes before turn it on again.
2. If the overload protector shuts off the motor frequently, you may have an unstable voltage situation.

Nebulizer Cleaning (JetNeb Kit)

All parts of the nebulizer, except tubing, should be cleaned according to the following instructions. Your physician and/or drug provider may specify a certain cleaning procedure. If so, follow their recommendations.

WARNING: Clean nebulizer after each use. Risk of bacteria growth, infection, illness from contamination due to dirty nebulizer.

WARNING!

Cleaning:

1. Ensure that the power cord is disconnected from the power outlet.
2. Disconnect nebulizer from tubing and compressor. Do not wash or submerge the air tubing.
3. Disassemble all parts of nebulizer by twisting the top and bottom, separating the pieces,
4. Clean nebulizer in warm soapy water.
5. Rinse all parts with clear water to remove any residue.
6. Allow the nebulizer and all parts to air dry completely and store in a clean plastic bag or container.

The 18260 JetNeb Nebulizer is intended for single patient use and should not be shared with other people. This unit should be replaced every thirty days.

CAUTION

Use and Cleaning of NEBKIT 500: refer to the user manual of NEBKIT 500 which is included in the accessory pack of NEBKIT 500 independently

Compressor Cleaning

1. With power switch in the "Off" position, unplug power cord from wall outlet.
2. Wipe outside of the compressor cabinet with a clean damp cloth to keep dust-free.

CAUTION: do not use any powdered cleaners or soap pads which may damage the finish.

CAUTION

WARNING: Do not submerge the compressor in water; doing so will result in compressor damage.

WARNING!

Filter Change


.....
DRIVE COMPRESSOR/NEBULIZER MODEL MQ5800, MQ5900, MQ5900P
.....

Change filter every 30 days or when filter turns gray.

CAUTION: Never operate the unit without a filter. Do not substitute cotton or any other materials for the filters. Do not attempt to clean the filters or reuse a filter.

To replace the filter, follow these steps:

1. Change the filter only when the unit is off, unplugged and completely cool.
2. Gently pull the filter cap out of the compressor.
3. Remove the filter from the filter cap.
4. Place a new filter in the filter cap.
5. Return the filter cap to its place, ensuring that it is firmly in place.

Maintenance

All maintenance must be performed by a qualified Drive provider or authorized service center.

Electric shock hazard. Do not remove compressor cabinet. All disassembly and maintenance must be done by a qualified Drive provider.

Size	7.5 in X 3.7 in X 5.7 in
Weight	3.0 lbs
Operating Voltage	120VAC, single phase, 60Hz
Current Consumption	2A (Max)
Storage/Transport Temperature Range	-13 to 158°F, (-25 to 70°C)
Storage/Transport Humidity	10 – 95%
Operating Humidity	10 – 95%
Operating Temperature Range	50 to 104°F, (10 to 40°C)
Compressor Flow Rate	≥9 L/minute
Operating Pressure Range	8~13 psi /55Kpa~90Kpa
Operating Noise Level	52 +\ -5 dB(A)
Service Life	1080 hours

18260 JetNeb Nebulizer

Min & Max filling volume	2~8 ml
Particle size (MMAD)	<5 μm
Nebulization rate	≥0.25 ml/min.
Disposable nebulizer	30 days

For specification of NEBKIT 500: refer to the user manual of NEBKIT 500, which is included in the accessory pack of NEBKIT 500 independently.

Equipment Classifications:

Protection against electric shock: Class II

Degree of protection against electric shock: Type BF

Equipment not suitable for use in the presence of flammable anaesthetic mixture with air.

Not intended for use with oxygen enriched environments such as oxygen hoods or tents.

*Conditions may vary based on altitude above sea level, barometric pressure, and temperature.

**Value determined with minimum flow rate through the nebulizer (20°C, sea level).

There are no hazardous materials used in the construction of this unit, therefore, there are no known risks associated with disposal of the equipment.

This product complies with the following electromagnetic compatibility standard:

IEC 60601-1-2

Please dispose of this nebulizer in accordance with local, municipal, state and federal law for electronic devices

Particle size (MMAD) – tested with Cascade impactor for 18260

Features	Drug	Result
Particle Size (MMAD) (µm)	Albuterol Sulfate (2.5 mg/3 ml)	3.40±0.44
	Ipratropium Bromide (0.5 mg/2.5 ml)	2.87±0.0.21
	Cromolyn Sodium (20 mg/2.0 ml)	3.43±0.81
Geometric Standard Deviation	Albuterol Sulfate (2.5 mg/3 ml)	3.07±0.35
	Ipratropium Bromide (0.5 mg/2.5 ml)	3.38±0.23
	Cromolyn Sodium (20 mg/2.0 ml)	3.74±0.97
Total Dose Delivered (µg)	Albuterol Sulfate (2.5 mg/3 ml)	759±62
	Ipratropium Bromide (0.5 mg/2.5 ml)	161±13
	Cromolyn Sodium (20 mg/2.0 ml)	3494±508
Total Respirable Dose (0.5-5 µm) (µg)	Albuterol Sulfate (2.5 mg/3 ml)	343±24
	Ipratropium Bromide (0.5 mg/2.5 ml)	72±9
	Cromolyn Sodium (20 mg/2.0 ml)	1528±169
Coarse Particle Dose (>4.7 µm) (µg)	Albuterol Sulfate (2.5 mg/3 ml)	390±49
	Ipratropium Bromide (0.5 mg/2.5 ml)	77±11
	Cromolyn Sodium (20 mg/2.0 ml)	1758±492
Fine Particle Dose (<4.7 µm) (µg)	Albuterol Sulfate (2.5 mg/3 ml)	369±18
	Ipratropium Bromide (0.5 mg/2.5 ml)	84±7
	Cromolyn Sodium (20 mg/2.0 ml)	1736±168
Ultra-Fine Particle Dose (<1.0 µm) (µg)	Albuterol Sulfate (2.5 mg/3 ml)	117±4
	Ipratropium Bromide (0.5 mg/2.5 ml)	29±4
	Cromolyn Sodium (20 mg/2.0 ml)	606±66

Recommended separation distances between portable and mobile RF communications equipment and the Compressor Nebulizer.

The Compressor Nebulizer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Compressor Nebulizer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Compressor Nebulizer as recommended below, according to the maximum output power of the communications equipment.

Recommended separation distances between portable and mobile RF communications equipment and the □

The Device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d = 1.2\sqrt{P}$	150 kHz to 80 MHz (in ISM and amateur radio bands) $d = 2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.7GHz $d = 2.3\sqrt{P}$
0,01	0.12	0.2	0.12	0.23
0,1	0.38	0.632	0.38	0.73
1	1.2	2	1.2	2.3
10	3.8	6.32	3.8	7.3
100	12	20	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic emissions

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Device is suitable for use in all establishments Including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class [B]	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic Immunity

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV input/output lines: ±1 kV	Power supply lines: ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV. line(s) to earth: ±2 kV. 100 kHz repetition frequency	line(s) to line(s): ±1 kV. 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5s	5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5s	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic Immunity

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance	
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150KHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Device ,including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d = 1.2\sqrt{P}$; $d = 2\sqrt{P}$	
Radiated RF IEC61000-4-3	10V/m	10V/m	80MHz to 800MHz: $d=1.2$ 800MHz to 2.7 GHz: $d = 2.3\sqrt{P}$	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the [XXXXX] is used exceeds the applicable RF compliance level above, the [XXXXX] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the [XXXXX].
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration - electromagnetic Immunity

The Device is intended for use in the electromagnetic environment specified below. The customer or the user of the Device should assure that it is used in such an environment

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 – 390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
	450	380 – 390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 – 787	LTE Band 13, 17	Pulse modulation b) 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
	870						
	930						
	1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
	5240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
	5240						
	5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E=6/d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

7 limited warranty

Five-Year Limited Warranty on Compressor Parts

Your Drive brand product is warranted to be free of defects in materials and workmanship for five years of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This 5 year Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.

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Drive Medical Canada Inc.
Donawa Lifescience Consulting Srl