SmartDose®

Mini Auto-Adjusting Oxygen Converter

drive™

EN CTOX-MN02

PRODUCT INFORMATION AND INSTRUCTIONS

- i** **NOTE:** Healthcare Equipment Providers shall provide the end user with instructions for use and instructions for determining gas cylinder contents.
- i** **NOTE:** The gauge face may be different on actual unit.

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TABLE OF CONTENTS

Important Information to Record.....	2	Checking for Leaks.....	15
Symbol Definitions.....	3	Operating Instructions.....	15
Important Safeguards, Dangers, Warnings and Cautions.....	3	Disassembly Instructions.....	17
Introduction.....	7	Caring for Your Conservar.....	17
Intended Use.....	8	Troubleshooting.....	18
Contraindications.....	8	Oxygen Cylinder Duration.....	20
Important Parts.....	9	Information for Home and Healthcare Providers.....	20
Standard Product.....	11	Disinfection Between Patients.....	20
Replacement Parts.....	11	Disinfection Intervals.....	22
Setting up Your Conservar.....	11	Maintenance.....	22
Assembly and Use.....	12	Calibration.....	23
Setup.....	13	Limited Warranty.....	23
Installing/Changing Alkaline Batteries.....	13	Disposal Instructions.....	23
Monitoring Battery Energy Level.....	13	Specifications.....	24
Assembly Instructions.....	14	Technical Description.....	25
Inspection Before Each Use.....	14	Electromagnetic Compatibility Information.....	26
		Pneumatic Diagram.....	27

IMPORTANT INFORMATION TO RECORD

Your Name: _____

Date You Received Your Unit: _____

Prescribed Oxygen Flow Setting:

- At Rest: _____

- During Exercise: _____

Home Care Provider's Name: _____

Home Care Provider's Phone Number: (_____) _____

Physician's Name: _____

Physician's Phone Number: (_____) _____

Notes: _____

SYMBOL DEFINITIONS

	Consult Instructions for Use		Humidity Limit
	Consult accompanying documents		Temperature Limits
	Prescription Only		Type BF Applied Parts
	No Smoking		Serial Number
	Use no oil or grease		Catalog Number
	Keep away from open flame		Manufacturer
	Use No Oil (on gauge face)		Date of Manufacture
cf	Continuous flow		Keep dry, protect the device from moisture
	MR Unsafe – Unsafe for Magnetic Resonance Environment		

DANGER!

Danger indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

WARNING!

Warning indicates a potentially hazardous situation which, if not avoided, will result in death or serious injury.

CAUTION!

Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

NOTE!

Gives useful tips, recommendations and information for efficient, trouble-free use.

NOTE: These symbols may be used in the user manual and also on the device, packaging and other material.

IMPORTANT SAFEGUARDS, DANGERS, WARNINGS AND CAUTIONS

Read and understand this manual before operating your oxygen conserving regulator. Important safeguards are indicated throughout this guide. Pay special attention to all safety information. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DANGER: Danger indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

- Oxygen will not burn; however, it does vigorously accelerate the burning of any flammable material.
-  Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula on bed coverings or chair cushions. If the oxygen conserver is turned on but not in use, the oxygen will vigorously accelerate the burning of any flammable material. Turn the oxygen conserver off when not in use.
-  Smoking during oxygen therapy is dangerous and is likely to result in burns or death. Do not allow smoking or open flames within the same room of the oxygen conserver or any oxygen carrying accessories.
- If you intend to smoke, you must always turn the oxygen conserver off, remove the cannula and leave the room where either the cannula or oxygen conserver is located.
-  There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.
- Avoid creation of any spark caused by any type of friction, near oxygen equipment.
-  Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
-  Never use oil, grease or petroleum-based products on or near the system to avoid the risk of fire and burns. Wash and dry hands before operating oxygen equipment. Use only water-based lotions or salves.
-  Do not lubricate replaceable fittings, connections, tubing or other accessories of the oxygen conserver to avoid the risk of fire and burns.
- Never use aerosol sprays near the equipment.

WARNING: Warning indicates a potentially hazardous situation which, if not avoided, will result in death or serious injury.

- Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using a CTOX-MN02 electronic oxygen conserver.
-  The use of oxygen requires a prescription from a physician or other appropriately licensed healthcare practitioner.
- The oxygen delivery setting has been determined for each patient individually with the configuration of the equipment to be used, including accessories.
- If you feel discomfort or are experiencing a medical emergency, seek medical assistance immediately to avoid harm.
- Geriatric or any other patient unable to communicate discomfort may require additional monitoring to avoid harm.
- Be sure to turn off the oxygen supply by pushing the CF flow control knob in.
- Oxygen flow is only present when flow indicator is clicked into place. No oxygen will flow if set between flow settings.
- Use of this device at an altitude above 9,800 ft. or outside temperature of 104° F is expected to adversely affect accurate delivery of oxygen therapy.
- Wind or strong drafts can adversely affect accurate delivery of oxygen therapy.
- This device is not intended to be used during sleep or by patients who breathe more than 40 breaths per minute or consistently fail to trigger equipment.
- The settings of this device may not correspond with continuous flow oxygen.

 **WARNING (Cont.)**

- The conserver must:
 - Be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
 - Be used with a specific combination of parts and accessories (i.e. cannulas, tubing) that meet the specifications of the conserver manufacturer, and that were used while your settings were determined.
- Your prescribed settings should be periodically reassessed by a physician or other appropriately licensed healthcare practitioner for effectiveness of therapy.
- The settings of other models or brands of oxygen therapy equipment do not correspond with the settings of the CTOX-MN02 Electronic Oxygen Conserving.
- Use only parts and accessories recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
- It is unsafe to use accessories, detachable parts and materials not described in the instructions for use, interconnect this equipment with other equipment not described in the instructions or modify the equipment.
- Cannula tubing can pose a strangulation risk and tripping hazard. Locate tubing away from areas of high foot traffic and keep children and pets away from cannula tubing.
- Always maintain a backup supply of oxygen (i.e. compressed cylinders and regulator).
- Not intended for use with any humidification devices.
- Not for use with an endotracheal tube or tracheostomy.
- Close supervision is necessary when the nasal cannula is used by impaired persons or near children.
- The proper placement of the nasal cannula in the nose is critical to the consistent operation of this equipment.
- Do not use cannula tubing that is longer than 7 ft. (2.13 m).

 **WARNING (Cont.)**

- Before using this unit, make certain that your hands are free of oil, grease and other contaminants, and that the cylinder is secure in an upright position.
- Lint and dust can contaminate the oxygen path and create a fire hazard. Keep all oxygen connections clean and free from lint and dust.
- Make sure the fill valve, if equipped, is clean and free of dirt, debris, oil or other contaminants.
- To reduce the risk of rapid oxygen recompression and fire, open the oxygen source valve slowly and completely so the pressure gauge moves slowly as it indicates the oxygen source pressure.
- Do not wet the area of the oxygen inlet connection, the disinfectant residue may create a risk of fire.
- Do not submerge the conserver in the disinfectant, wipe disinfect only.
- Do not attempt to open or remove enclosure. There are no user-serviceable internal components.
- MR Unsafe. Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the conserving device or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.
- Choking Hazard - The 870 post seal washer may pose a choking hazard. Small parts not for children under 3 years or any individuals who have a tendency to place inedible objects in their mouths.

WARNING (Cont.)

- Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the conserving device. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, discontinue use. Contact your home care provider. This device is suitable for use in home and healthcare environments except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of Electromagnetic DISTURBANCES is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the conserving device. Otherwise, degradation of the performance of this equipment could result.

CAUTION: Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

- Prevent water or any other liquid substances from entering the unit.
- Prevent dust, lint or any small particles from entering the unit, as these may degrade performance or cause oxygen leaks.
- Do not expose the unit to extreme temperatures or heat sources such as fireplaces and radiant heaters.
- Do not use humidifier bottles.
- Do not use if leaking or damaged.
- Refer repairs to authorized service personnel.
- Do not use the flow outlet for driving any medical equipment.
- Oxygen supplied by this equipment is supplemental only and is not intended for life support applications. The CTOX-MN02 electronic oxygen conserver should not be used to supply anything other than medical oxygen.
- Remember that in continuous flow mode, the oxygen will be consumed at a much faster rate. Switch to another oxygen source before depleting the oxygen cylinder. The continuous flow function on the CTOX-MN02 electronic oxygen conserver is factory preset at 2 LPM.
- Checking for leaks should be performed before each use.
- To help prevent possible damage to the unit and to maintain its cleanliness, keep the CTOX-MN02 electronic oxygen conserver in a carrying bag. Several bags are available for use with different cylinder sizes and configurations: contact your Home Care provider for more information on carry bags.
- Combinations with medical devices other than recommended can alter the performance of the equipment (example: combinations with reservoir cannula, extended length tubing).

SAVE THESE INSTRUCTIONS

INTRODUCTION

This manual provides information necessary to operate the CTOX-MN02 electronic oxygen conserver. This conserver can be used with CGA 870 post-valve cylinder [see Fig. A] at home or away from home to provide your specific oxygen requirements. It requires two (2) 1.5 volt "AA" alkaline batteries for operation.

PURPOSE:

The SmartDose CTOX-MN02 electronic oxygen conserver includes a combination of a low-pressure regulator and an oxygen conserver. It is designed for use with a cylinder as an ambulatory oxygen system and is capable of delivering a specific amount of supplemental oxygen at the optimal point in the breathing cycle. The SmartDose increases the efficiency in the delivery of oxygen, maximizing the beneficial effects and eliminating unnecessary oxygen waste.

FUNCTION:

When we breathe, approximately one-third of the time is spent inhaling and two-thirds exhaling. As a result, oxygen delivered by continuous flow is wasted during exhalation. By eliminating oxygen flow during exhalation, a two-thirds saving is possible. Additionally, the oxygen available during the very first part of inhalation contributes most to meeting oxygen needs. The CTOX-MN02 electronic oxygen conserver takes advantage of these facts to provide maximum efficiency in the delivery of oxygen. This device is designed to be an integral component of a lightweight, long-lasting supplemental ambulatory oxygen system.

USER QUALIFICATION:

Prior to beginning therapy, patients must be given instruction by qualified personnel in how to operate this unit.

Functions that are frequently used include:

- Installing the device by connecting to a portable oxygen cylinder post valve.
- Selecting prescribed oxygen setting on the selector dial.
- Attaching standard oxygen cannula and tubing to gas port and securing cannula to user.
- Monitoring the pressure gauge for remaining oxygen in cylinder.
- Disconnecting the device from the portable oxygen cylinder.


FIGURE A
Post-Valve Cylinder

INTENDED USE:

The SmartDose Oxygen Conserver is intended as a delivery device for medical-grade oxygen from high-pressure oxygen cylinders. This is an ambulatory device, which allows patients to ambulate longer than they would with a continuous flow regulator on the same cylinder.

The SmartDose Oxygen Conserver is intended to be used in the hospital, healthcare facilities, or home care environments.

The patient is the intended operator.

 **CAUTION:** Oxygen supplied by this equipment is supplemental only and is not intended for life support applications. The CTOX-MN02 electronic oxygen conserver should not be used to supply anything other than medical oxygen.

CONTRAINDICATIONS:

- This device is not to be used for life support applications.
- This device is not intended for use by patients who breathe more than 40 breaths per minute.
- This device is not to be used by patients who consistently fail to trigger the device.
- The device is not for use while asleep.

 **NOTE:** The CTOX-MN02 electronic oxygen conserver will be referred to as the SmartDose in the remainder of the guide.

IMPORTANT PARTS

- 1. Cylinder Attachment Knob:** Used to attach the Converter to the cylinder.
- 2. SmartDose Information Signal:** One green light illuminates one flash on each breath when in SmartDose mode to show the current setting.
- 3. SmartDose On/Adjust Button:** Use this to power on and adjust the patient setting. Push and release to turn on. Push and hold to change settings.
- 4. Low Battery Information Signal:** Flashes red when batteries are low and should be replaced as soon as possible. Solid red means batteries are dead.
- 5. Cannula Fitting:** Use this fitting to attach the cannula tubing.
- 6. Oxygen Contents Gauge:** Indicates the remaining pressure in the oxygen cylinder. When this gauge falls into the red section, you should switch to a new cylinder.
- 7. Continuous Flow (CF) Knob:** Pull knob out to switch unit to 2 LPM CF oxygen, does not require batteries, must be in to initiate SmartDose.
- 8. Battery Compartment:** Battery holder for 2 "AA" alkaline batteries. Note polarity when inserting new batteries.


FIGURE B

The CF mode is a mechanical override designed as a safety feature in case of battery power failure. When in the CF mode, the user can still change the LED display lights by pressing the Selector Button, but it has no effect on the oxygen output.

CAUTION: Remember that in continuous flow mode, the oxygen will be consumed at a much faster rate. Switch to another oxygen source before depleting the oxygen cylinder. The continuous flow function on the SmartDose converter is factory preset at 2 LPM.

The normal operating mode is with the CF knob pushed in and pulses oxygen at the rate set by the Selector Button.

NOTE: When using the SmartDose in Continuous Flow mode, the information signal light will NOT light up green or red with each breath. Turning the On/Adjust button all the way down will not completely stop the flow of oxygen. To stop the flow of oxygen push the CF switch in & turn the cylinder valve off.

CAUTION: In the event that it is necessary to operate the unit in the continuous flow mode, DO NOT obstruct the flow of oxygen from the Oxygen Supply Outlet by placing your finger over the outlet or blocking the flow through the oxygen tubing in any way. Doing so may render the unit inoperable and/or damage the sensor in the unit.

IMPORTANT PARTS (Cont.)


WARNING: Fire Hazard

Lint and dust can contaminate the oxygen path and create a fire hazard. Keep all oxygen connections clean and free from lint and dust.

SMARTDOSE MINI AUTO-ADJUSTING OXYGEN CONSERVER DELIVERY MODES

The SmartDose Mini Auto-Adjusting Oxygen Conserver has two modes: SmartDose and Continuous Flow. To receive the benefits of SmartDose technology, set your SmartDose Mini Auto-Adjusting Oxygen Conserver to SmartDose mode. Continuous flow mode is offered as a backup option in case your batteries run out of power. Continuous Flow mode provides a steady flow of oxygen until you can replace the batteries; it does not offer the benefits of SmartDose technology. Also, be aware that your oxygen supply will be depleted rapidly when set to operate in Continuous Flow mode.

**SmartDose mode
CF knob pushed IN**

**SmartDose mode CF
knob pushed OUT**

With the conserver in SmartDose mode, SmartDose technology detects changes in your breath rate and automatically adjusts the oxygen dose using a unique valving system. In this way, the Conserver automatically adapts the oxygen flow to meet your needs as your activity level goes up and down throughout the day. After you have set the flow rate to your prescribed dose, there is no need to manually change the oxygen flow rate during the day— SmartDose technology does it for you.

WARNING: Use only a manufacturer-specified seal washer. Other seal washers may not be oxygen compatible and may cause an oxygen leak, creating an increased fire risk.

STANDARD PRODUCT

ORDER NUMBER	DESCRIPTION
CTOX-MN02	Electronic oxygen conserving device with CGA 870 post connection

REPLACEMENT PARTS

ORDER NUMBER	DESCRIPTION
SEAL-100	870 Post Seal Washer (Gasket)

SETTING UP YOUR CONSERVER

WARNING: Fire Hazard

Before using this unit, make certain that your hands are free of oil, grease and other contaminants, and that the cylinder is secure in an upright position.

NON-PORTABLE USE:

The SmartDose is designed to extend the life of portable oxygen supplies when away from the primary source. While the conserver may be used with stationary oxygen sources, the unit should be used only while awake. The conserver is not intended for use while asleep.

Conditions in the home healthcare environment that can affect the performance of your device:

WARNING: Fire Hazard

- Lint and dust can contaminate the oxygen path and create a fire hazard. Keep all oxygen connections clean and free from lint and dust.
- Degraded sensors due to rough handling, extreme temperatures or humidity can affect the device's ability to dose oxygen. To avoid this, operate the device within the rated environmental conditions listed in the Specifications section of this guide.

ASSEMBLY AND USE

WARNING: Fire Hazard

- Make certain that your hands are free of oil, grease and other contaminants.
- Inspect the unit to ensure the seal washer (gasket) is in good working condition and attached to the inlet nozzle.
- Secure the cylinder in an upright position.

WARNING: Fire Hazard

- Inspect the valve of the cylinder and the conserver to ensure they are free of contaminants and damage. If any indication of damage or contamination is detected, **DO NOT** use the equipment and contact your Home Care Provider.

- **WARNING:** Use **ONLY** a manufacturer-specified seal washer (gasket). An incorrect seal washer (gasket) may not be oxygen compatible or may cause an oxygen leak, creating an increased fire risk. **DO NOT** use the device if the manufacturer-specified seal washer (gasket) is missing.

SETUP

INSTALLING/CHANGING ALKALINE BATTERIES:

The SmartDose uses advanced technology that prolongs battery life through efficient power use. With normal use of four (4) hours per day at any setting, your batteries should last a minimum of 1 year.

Two (2) "AA" alkaline batteries are included in your SmartDose. Follow these easy steps to install or change batteries when needed:

STEP 1: Remove the battery cover by pressing the clip on the top of the unit and pull the cover out.

STEP 2: Remove old batteries from the compartment and discard properly.

STEP 3: Insert 2 "AA" Alkaline with the polarity as indicated on the holder.

STEP 4: Insert the battery cover back onto the unit and push until the retaining clip snaps into place.


FIGURE C

i NOTE: The SmartDose is packaged with batteries inside. A special seal is used at the positive terminal to prevent battery oxidation. If it has not already been removed by your Home Care Provider, be sure to remove the seal before using the unit for the first time.

MONITORING BATTERY ENERGY LEVEL:

The SmartDose is equipped with a visual low battery information signal that alerts you when a battery change is needed. The following signals will occur as your batteries become depleted:

- Blinking red light: Low Battery energy level. Replace batteries.
- Solid red light: Batteries are dead.

Disposal of Batteries: Do not dispose of used batteries in the household waste. Contact a public waste disposal authority for proper disposal instructions.

ASSEMBLY INSTRUCTIONS:

STEP 1: Loosen the cylinder attachment knob.

STEP 2: Lower the SmartDose conserver over any CGA-870 post-valve cylinder with the alignment pins toward the holes on the cylinder neck [see Fig. D].

STEP 3: Line up the two pins and the seal washer (gasket) with the corresponding holes on the cylinder post valve.

STEP 4: While holding the unit in place, tighten the cylinder attachment knob by turning clockwise [see Fig. B].

i NOTE: The cylinder attachment knob should be aligned with the indentation on the post valve.

i NOTE: Tighten only by hand. The use of a tool to tighten the knob may damage the unit.

STEP 5: Attach a standard cannula (7 ft. [2.13 m] or less in length) to the oxygen supply outlet. See Fig. E for an illustration demonstrating the proper positioning of the nasal cannula.

INSPECTION BEFORE EACH USE:

STEP 1: Examine the outside of the cylinder and regulator assembly for debris, oils or grease and noticeable signs of damage, dents or dings, corrosion, excessive heat or fire damage, pressure leak or any other sign of damage that might cause a cylinder to be unacceptable or unsafe for use.

! WARNING: Fire Hazard
If any of these conditions exist, DO NOT use the cylinder and contact your healthcare provider for a replacement cylinder.

STEP 2: Check the continuous flow (CF) knob to be sure it moves in and out freely.


FIGURE D

CHECKING FOR LEAKS:

CAUTION: This should be performed before each use.

STEP 1: Close the valve on the oxygen cylinder by turning the valve clockwise.

STEP 2: Depressurize the conserver by inhaling several times using the nasal cannula. The gauge indicator should drop to zero.

STEP 3: Check that all tube connections are tight. If necessary, tighten them **by hand**.

CAUTION: Do not use a tool.

STEP 4: Ensure that the unit is set to the "SmartDose" position (CF knob pushed in).

STEP 5: Slowly open the valve on the oxygen cylinder by turning the valve counter clockwise until the needle in the gauge indicator no longer moves.

STEP 6: Close the oxygen cylinder valve again by turning the valve clockwise.

STEP 7: Observe the needle in the gauge indicator for approximately one minute.

- If the needle remains in its position, there is no leak and the device is ready for use.
- However, if the contents indicator shows a continuous decrease in pressure, there is a leak in the system. In this event, contact your Home Care Provider.

OPERATING INSTRUCTIONS:

STEP 1: Make sure that the SmartDose CF knob is set to the "SmartDose" position (CF knob pushed in) before opening the oxygen source valve.

STEP 2: Slowly open the valve on the oxygen cylinder by turning the valve counter-clockwise until the needle in the gauge indicator no longer moves.

WARNING: Fire Hazard

To reduce the risk of rapid oxygen recompression and fire, open the oxygen source valve slowly and completely so the pressure gauge moves slowly as it indicates the source pressure.

STEP 3: If a leak is present (see "Checking for Leaks" section above), close the oxygen cylinder valve, check the seal washer and reinstall. If the leak persists, **DO NOT USE THE EQUIPMENT**. Contact your Home Care Provider for repair or replacement.

STEP 4: Check the oxygen pressure gauge to verify that the oxygen cylinder pressure is within the operating range.

STEP 5: Turn on the SmartDose by pressing and holding the On/Adjust button until all of the numerical LEDs (1 through 5) light up. The SmartDose indicates it is ready by pulsing three times.

i NOTE: To check your current setting, press and release the SmartDose On/Adjust Button. Set the information signal light to your prescribed setting (1 to 5). Be sure that the CF (Continuous Flow/SmartDose switch is pushed in as shown, confirming that your SmartDose conserver is in SmartDose mode (See Fig. B).

STEP 6: Ensure cannula is properly attached to the SmartDose conserver.

STEP 7: Properly position the cannula with the nasal prongs facing upward and insert the prongs into nose. Wrap the cannula tubing over the ears and position in front of body(see Fig. E).

FIGURE E

The conserver will now start to deliver oxygen. The amount of oxygen delivered per pulse is determined by the setting. A sound may be heard each time the unit delivers a pulse of oxygen. Adequate saturation will be achieved because of the precise time in the breathing cycle in which the pulse of oxygen is delivered.

! CAUTION: To help prevent possible damage to the unit and to maintain its cleanliness, keep the SmartDose in a carrying bag. Several bags are available for use with different cylinder sizes and configurations: contact your Home Care provider for more information about carry bags.

i NOTE: If the SmartDose is in Continuous Flow mode, 2 LPM of oxygen will still be delivered. The conserver must be powered off & the CF switch pushed in to completely turn off the flow of oxygen.

STEP 8: When finished using the system, close the oxygen cylinder valve and continue breathing through the nasal cannula until no further oxygen is detected and the gauge reads "empty" or "zero (0)".

STEP 9: Remove the nasal cannula.

i NOTE: When you are finished using the SmartDose simply remove the cannula tubing. After 5 minutes of inactivity, the unit will power itself down. The unit quickly opens and closes the valves 3 times thirty seconds before powering down to alert you that it is about to shut down. In the event that you have simply been mouth breathing, this should alert you that the unit has not seen a breath and will be shutting down. Close cylinder valve when not in use.

STEP 10: When not in use, store in a clean, dry location.

DISASSEMBLY INSTRUCTIONS:

To remove the SmartDose from the oxygen cylinder:

STEP 1: Close the oxygen cylinder valve by turning the valve clockwise.

STEP 2: Deplete residual pressure in the conserver by continuing to breathe through the nasal cannula or by setting the conserver to continuous flow mode until no further oxygen is detected and the gauge reads "empty" or "zero (0)".

STEP 3: Remove the nasal cannula from the SmartDose and push the CF knob in to SmartDose mode.

STEP 4: While holding the unit in place, slowly loosen the cylinder attachment knob by turning counter-clockwise (see Fig. B).

STEP 5: Pull the unit away from the cylinder post valve until the two alignment pins are out of the holes on the cylinder post valve.

STEP 6: Lift the unit off of the cylinder post valve. Use caution while lifting the unit off of the post valve, as the alignment pins can damage sealing surfaces on the post, increasing the chance of leakage.

STEP 7: Store in a clean, dry location.

CARING FOR YOUR CONSERVER

The SmartDose should be kept clean and free from moisture and dust, as well as extreme temperature. Clean the SmartDose once a week by wiping the outer surface and oxygen outlet port with a clean, dry, lint-free cloth. Do not use oil or products containing oil.

Do not expose the unit to water, such as when bathing or swimming. It is advisable to keep the system in its carrying bag to afford a degree of protection. Clean the outside of the unit periodically with a clean, lint-free cloth. To remove heavy soil, wipe unit using a lint-free cloth, dampened with clean water, then dry unit with clean, dry, lint-free cloth. Pay special attention to the oxygen outlet to make sure it remains free of dust, etc.

WARNING: There are no user serviceable parts. Do not attempt to service the device while it is in use.

TROUBLESHOOTING

PROBLEM	PROBABLE CAUSE	SOLUTION
Oxygen is not being delivered, even though green light is flashing on the Conserver each time I inhale.	Oxygen cylinder is empty.	Check contents gauge on unit. If pressure is in the red, switch cylinders.
	Oxygen supply is not turned on.	Open the oxygen flow by opening the valve on the top of the cylinder.
Use times are different from those stated in the literature.	SmartDose responds to your breath rate and volumes will be increased as your need increases. Time shown are estimates and will vary from patient to patient and day to day.	SmartDose is working properly.
	Leak in system.	Check connection to cylinder. If leak is apparent, contact your oxygen provider for new seal.
SmartDose will not pulse - no green light when I inhale.	Cannula is not on properly and the unit cannot sense your inhale.	Check the cannula tubing and ensure that it is comfortably in your nose.
	The unit is not turned on or has powered down.	Press and hold Power/Adjust or On/Adjust button to turn the unit on.
	Batteries are depleted or not installed.	Install 2 x "AA" new Alkaline batteries.
	Mouth breathing.	Inhale through nose.
	Unit is in Continuous Flow (CF) mode.	Push CF knob in all the way to switch to SmartDose mode.

TROUBLESHOOTING (Cont.)

PROBLEM	PROBABLE CAUSE	SOLUTION
SmartDose works fine for some time and then sensitivity seems to drift or unit stops working.	Cannula has too much back pressure (pediatric).	Replace cannula with standard cannula.
	Pressure sensor is not zeroed.	Power off unit by removing batteries. Remove cannula from nose. Reinstall batteries, power on unit, insert cannula and resume normal breathing.
Short battery life.	Non-alkaline batteries are used.	Make sure the batteries inside the unit are alkaline.
	Batteries are faulty.	Replace batteries. (See "Installing Alkaline Batteries" section).

Non-functioning units are subject to warranty provisions and the manufacturer repair/return policy. If necessary, call your Home Care Provider.

 **NOTE:** Do not attempt to open the unit. If the unit is opened or tampered with, the warranty is void.

OXYGEN CYLINDER DURATION

Because the total delivery of oxygen via the SmartDose is related to breathing rates, it is user adaptive. The total oxygen delivered per minute will automatically adjust with user need. For example, at all settings, twice as much oxygen per minute will be delivered if one breathes twenty (20) times per minute as compared with ten (10) times per minute. Please refer to the table below as a guide.

SmartDose CTOX-MN02	SETTING	1	2	3	4	5	CF 2 LPM
Cylinder Type	Cylinder Volume (Liters)	Estimated Cylinder Duration in Hours (based on 20 breaths per minute)					
M6(B)	164 liters	8.3	4.1	2.7	2.0	1.6	1.4
ML6	171 liters	8.6	4.3	2.8	2.1	1.7	1.4
M9(C)	246 liters	12.4	6.1	4.1	3.0	2.4	2.1
D	425 liters	21.5	10.6	7.0	5.2	4.2	3.5
E	680 liters	34.3	16.9	11.2	8.4	6.7	5.7

INFORMATION FOR HOME AND HEALTHCARE PROVIDERS

DISINFECTION BETWEEN PATIENTS:

Drive Medical recommends that at least the following procedures be completed by an appropriately trained individual (i.e., biomedical technician, equipment technician, respiratory therapist, nurse) between uses by different patients.

STEP 1: Dispose of nasal cannula.

STEP 2: Cleaning Instructions:

- Clean the exterior of the device with:
 - a clean, lint-free cloth to remove surface dirt and soil.
 - a clean, lint-free cloth dampened with water to remove heavy soil.
 - a soft-bristled brush dampened with water to remove stubborn soil.
- Dry the device with a clean, lint-free cloth.

 **NOTE:** Ensure that the oxygen inlet remains free of dust, etc.

STEP 3: Disinfection Instructions:

- Saturate a clean, lint-free cloth with a suitable disinfecting agent (e.g., 1:5 chlorine bleach [5.25%] and water solution) and wipe the exterior of the device.

NOTE: Ensure that the device remains visibly wet for ten minutes, with the solution reaching all crevices and hard-to-reach areas.

WARNING: Do not wet the area of the oxygen inlet connection, the disinfectant residue may create a risk of fire. Do not submerge the conserver in the disinfectant, wipe disinfect only.

STEP 4: Allow all parts to dry completely before checking for proper operation.

STEP 5: Check the oxygen conserver for proper operation:

- Press and hold the Selector Button and ensure that each pulse setting can be selected and that each Dosage LED illuminates.
- Attach device.
- Ensure pressure gauge is working.
- Ensure CF flow is working.
- Set device to each setting, and verify that oxygen pulses.

If the device is operating properly, it is ready to be used by another patient. If the device is not operating properly, the provider should return the device to Drive Medical for service.

NOTE: There is no portion of the gas pathways through the conserver that could be contaminated with body fluids or expired gas under normal or single fault conditions.

DISINFECTION INTERVALS

 **WARNING:** Do not attempt to open or remove enclosure. There are no user-serviceable internal components.

 **NOTE:** The disinfection process can only be performed by an appropriately trained individual.

PART/ACCESSORY	RECOMMENDED DISINFECTION INTERVAL	COMPATIBLE DISINFECTION METHOD
Outside surface of the conserver and controls	Between patients	1:5 chlorine bleach (5.25%) and water solution
Oxygen Tubing, nasal cannula	Do not clean, replace between patients	N/A

MAINTENANCE

The oxygen conserver does not need maintenance through the expected service life of the product. The device has an expected service life of 5 years, which is based on 1460 hours/year of use under normal operating and storage conditions.

PREVENTIVE MAINTENANCE

Before attaching the conserver to an oxygen cylinder:

- Examine the seal washer and replace if it is damaged or contaminated with dirt, debris, oil or other contaminants.
- Inspect the threads on the cylinder attachment knob, make sure that the knob turns freely through its full range of adjustment.

To replace the SEAL-100 Seal Washer:

1. Remove the old seal washer by lifting the edge of the brass washer and removing the seal washer from the oxygen inlet connection, discard old seal washer.
2. Wash and dry hands before handling the new seal washer.
3. Place the new seal washer over the oxygen inlet connection and press it against the body of the device to seat the seal washer in place.

CALIBRATION

No calibration is necessary through the expected life of the product.

LIMITED WARRANTY

The SmartDose conserver has been carefully manufactured and inspected and is warranted to be free from defects in workmanship and materials. Under this warranty, Drive Medical's obligation shall be limited to the replacement or repair of any such units or parts that prove, by Drive Medical's inspection, to be defective within two years from the date of purchase. Any abuse, operation other than the intended use of the product, negligence, accident or repair by other than authorized service professionals shall immediately void this warranty. This warranty does not extend to spare parts or accessories.

Drive Medical will not accept damages or charges for labor, parts or expenses incurred in making field repairs, except upon written authorization prior to such action.

The foregoing warranty is exclusive and in lieu of all other express warranties. Implied warranties, if any, including but not limited to the implied warranties of merchantability and fitness for a particular purpose, shall not extend beyond the duration of the express warranty provided herein. In no event shall Drive Medical be liable for loss of use or profit or other collateral, special or consequential damages.

DISPOSAL INSTRUCTIONS

The equipment, including accessories and internal components, does not belong in your regular household waste. Such equipment is manufactured from high-grade materials and can be recycled and reused. The plastic components should be recycled as plastics.

The device should be destroyed or made unusable prior to disposal, example: remove the internal PC board.

Dispose in accordance with local ordinances.

Disposal of Batteries: Do not dispose of used batteries in the household waste. Contact a public waste disposal authority for proper disposal instructions.

SPECIFICATIONS

DELIVERED OXYGEN VOLUME PER BREATH (ml)						
Setting	Breath Rate (BPM)					
	15	20	25	30	35	40
1	15	15	15	15	15	16
2	30	30	27	27	27	28
3	44	39	39	40	40	40
4	52	52	53	52	50	43
5	65	65	65	60	50	43

Deviation of delivered oxygen per breath =
 +/- (1 ml/breath + 10% of the set initial value
 per breath)
 Tested per ISO 80601-2-67:2014 Sec 201.12.1.101

Continuous Flow Rate Setting cf =
 2 LPM +/- 0.4 LPM

Tolerance on measurements = +/- 15%
 (including measurement uncertainty)

All volumes and flows are expressed at STPD
 (standard temperature and pressure, dry)

Continuous Flow Emergency Bypass System Setting:	Factory preset at 2 LPM
Regulator:	Brass high-pressure with aluminum low-pressure materials
Dimensions (LxHxW) :	Approximately 12.7 cm (5.0 in) x 12.7 cm (5.0 in) x 9.5 cm (3.8 in) conserver only
Weight:	Approximately 17.6 ounces (500 grams) conserver with batteries
Operating Temperature:	41°F to 104°F (5°C to 40°C)
Operating Relative Humidity:	15% to 93% non-condensing
Operating Altitude:	-1,200 to 9,800 ft.
Operating Atmospheric Pressure:	1,060 to 700 hPa
Operating Inlet Pressure:	500 to 3000 psig (34.5 to 206.8 bar)
Storage/Transportation:	-13°F to 158°F (-25°C to 70°C) RH up to 93% non-condensing
Shock:	IEC 60601-1-11:2010
Vibration:	IEC 60601-1-11:2010
Pneumatic Inspiratory Trigger Sensitivity:	0.03-0.56 cm H ₂ O sensitivity: (0.01-0.22 in H ₂ O) (under test conditions in ISO 80601-2-67:2014 sec. 201.12.1.101)
Breathing Frequency Range:	5 to 40 bpm
Ingress Protection Rating:	Keep Dry
Batteries:	(2) x 1.5 volt Alkaline type AA
Low Battery Information Signal-Blinking LED:	Low Battery Energy Level, replace batteries

Technical Description

How the device functions:

The SmartDose includes a combination of a low-pressure regulator and an oxygen conserver. It is designed for use with a cylinder as an ambulatory oxygen system and is capable of delivering a precise amount of supplemental oxygen at the optimal point in the breathing cycle. The SmartDose greatly increases the efficiency in the delivery of oxygen, maximizing the beneficial effects and eliminating unnecessary oxygen waste. The SmartDose is designed to sense the precise moment in the patient's breathing cycle when it has to deliver the oxygen pulse. It accomplishes its intended purpose by providing supplemental oxygen to a patient requiring nasal oxygen therapy in a manner that yields a high level of efficiency in the consumption of the oxygen administered. The SmartDose conserver's responsive and rapid oxygen delivery time increases the oxygen concentration provided to the patient, thus increasing the ability to better maintain patient saturation levels.

Principles of dosing, timing, triggering and settings.

The normal volume of oxygen that SmartDose delivers is 16 ml per setting, so that a setting of 2 delivers 32 ml of oxygen on each breath, in the first 2/3 of the inhale cycle. SmartDose monitors a patient's breathing patterns and automatically turns up the oxygen volume by one setting during exertion, and if breath rate continues to rise, by another setting. Many patients turn their settings up and down throughout the day, but this self administered approach is inconsistent and by the time most patients feel the need to turn it up, their oxygen levels have already dropped. SmartDose responds more quickly and consistently. SmartDose recognizes when exertion stops and lowers the volume back down to normal levels thus saving oxygen. Again, this is often an improvement over manual adjustments since patients often forget to turn the volume down after exerting. Most patients are able to be titrated at lower base settings due to SmartDose's efficient delivery and auto-adjusting features. This extends the use time of your oxygen supply. SmartDose also uses a multi-valve system that has multiple delivery flow rates. The unit can deliver a volume at two flow choices. The unit constantly monitors the patient's breath rate and determines the lowest flow option to "fit" the oxygen into the first 2/3 of the inhalation cycle. By delivering at lower flows, the unit is more comfortable, and quieter than units that deliver at higher flows all the time.

 **NOTE:** The components in the medical gas path and the nasal cannula are considered applied parts.

 **CAUTION:** Combinations with medical devices other than recommended can alter the performance of the equipment (example: combinations with reservoir cannula, extended length tubing).

Electromagnetic Compatibility Information

Test Description	Specification	Notes	Results
CISPR11 & IEC 60601-1-2:2014 (4th Edition)			
Radiated Emissions	EN 55011:2009 + A1:2010 FCC Part 15, Subpart B	Class B 30MHz – 1GHz	Complies
Electrostatic Discharge Immunity	EN 61000-4-2:2008	±15kV (Air) 8kV (Contact)	Complies
Radiated Electromagnetic Field Immunity	EN 61000-4-3:2010	80MHz to 2.7GHz @ 10V/m; 80% AM at 1kHz	Complies
Proximity Fields from RF Wireless	EN 61000-4-3:2010	385Mhz – 5785Mhz	Complies
Magnetic Field Immunity	EN 61000-4-8:2009	30A/m @ 50Hz & 60Hz	Complies

FCC Statement


i NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced technician for help.

Pneumatic Diagram

Pneumatic Diagram depicts device in conserving selection.

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Inovo, Inc.
DeVilbiss Healthcare Pty Limited
Printed in the U.S.A