

## Non-Tilt Pneumatic Cylinder Overbed Table H

HCPCS:



FIGURE A



- Pneumatic cylinder lift mechanism infinite height adjustment from 27" 42"
- ¾" x 18" x 32" high pressure laminate top with T- molding
- 2" casters for easy maneuvering
- Column can be mounted on left or right side
- Available with "H" (Figure A) or "U" base in chrome or opal powder coat
- Optional dual mirror vanity drawer
- Base and columns available in additional powder coat finishes
- Available on drop ship basis only

| PRODUCTS    | Description                      |
|-------------|----------------------------------|
| 16008-HCBO  | OBT-PneumCyl-BO-H-Chrm-SpclOrdr  |
| 16008-HCBOV | OBT-PnmCyl-BO-H-Ch-Van-SpclOrdr  |
| 16008-HCDM  | OBT-PneumCyl-FM-H-Chrm-SpclOrdr  |
| 16008-HCDMV | OBT-PnmCyl-FM-H-Ch-Van-SpclOrdr  |
| 16008-HCEO  | OBT-PneumCyl-EO-H-Chrm-SpclOrdr  |
| 16008-HCEOV | OBT-PnmCyl-EO-H-Ch-Van-SpclOrdr  |
| 16008-HCFM  | OBT-PneumCyl-FM-H-Chrm-SpclOrdr  |
| 16008-HCFMV | OBT-PnmCyl-FM-H-Ch-Van-SpclOrdr  |
| 16008-HCGW  | OBT-PneumCyl-GW-H-Chrm-SpclOrdr  |
| 16008-HCGWV | OBT-PnmCyl-GW-H-Ch-Van-SpclOrdr  |
| 16008-HCNT  | OBT-Pneumcyl-NT-H-Chrm-SpclOrdr  |
| 16008-HCNTV | OBT-PnmCyl-NT-H-Ch-Van-SpclOrdr  |
| 16008-HCQO  | OBT-PneumCyl-QO-H-Chrm-SpclOrdr  |
| 16008-HCQOV | OBT-PnmCyl-QO-H-Ch-Van-SpclOrdr  |
| 16008-HCSO  | OBT-PneumCyl-SO-H-Chrm-SpclOrder |
| 16008-HCSOV | OBT-PnmCyl-SO-H-Ch-Van-SpclOrder |
| 16008-HCWC  | OBT-PneumCyl-WC-H-Chrm-SpclOrdr  |
| 16008-HCWCV | OBT-PnmCyl-WC-H-Ch-Van-SpclOrdr  |
| 16008-HOBO  | OBT-PnmCyl-BO-H-Op-SpclOrdr      |
| 16008-HOBOV | OBT-PnmCyl-BO-H-Op-Van-SpclOrdr  |
| 16008-HODM  | OBT-PnmCyl-FM-H-Op-SpclOrdr      |
| 16008-HODMV | OBT-PnmCyl-FM-H-Op-Van-SpclOrdr  |
| 16008-HOEO  | OBT-PnmCyl-EO-H-Op-SpclOrdr      |
| 16008-HOEOV | OBT-PnmCyl-EO-H-Op-Van-SpclOrdr  |
| 16008-HOFM  | OBT-PnmCyl-FM-H-Op-SpclOrd       |

## PRODUCTS

## Description

|             |                                  |
|-------------|----------------------------------|
| 16008-HOFMV | OBT-PnmCyl-FM-H-Op-Van-SpclOrdr  |
| 16008-HOGW  | OBT-PnmCyl-GW-H-Opal             |
| 16008-HOGWV | OBT-PnmCyl-GW-H-Opal-Vanity      |
| 16008-HONT  | OBT-PnmCyl-NT-H-Op-SpclOrdr      |
| 16008-HONTV | OBT-PnmCyl-NT-H-Op-Van-SpclOrdr  |
| 16008-HOQO  | OBT-PnmCyl-QO-H-Op-SpclOrdr      |
| 16008-HOQOV | OBT-PnmCyl-QO-H-Op-Van-SpclOrdr  |
| 16008-HOSO  | OBT-PnmCyl-SO-H-Op-SpclOrder     |
| 16008-HOSOV | OBT-PnmCyl-SO-H-Op-Van-SpclOrder |
| 16008-HOWC  | OBT-PnmCyl-WC-H-Op-SpclOrdr      |
| 16008-HOWCV | OBT-PnmCyl-WC-H-Op-Van-SpclOrdr  |