



HCPCS:











FIGURE A

FIGURE B

FIGURE C

FIGURE D



FIGURE E

- Height adjustable in 1" increments
- Accessories mount tool free and can be adjusted according to the user's needs
- Can be used in anterior or posterior position
- Medium and Large frames fold for easier transportation
- Swivel casters can be locked in non-swivel position
- Directional lever allows for both forward and reverse mobility, or forward mobility only
- Wheel locks on all 4 casters to keep gait trainer stationary during transfers
- Small and Medium sizes available in blue and green
- Large size available in grey

PRODUCTS	Description	
GT1000-2GB	Moxie GT Gait Trainer, Small, Blue	
GT2000-2GB	Moxie GT Gait Trainer, Medium, Blue	
GT2000-2GG	Moxie GT Gait Trainer, Medium, Green	
GT3000-2GGY	Moxie GT Gait Trainer, Large, Grey	
SPECIFICATIONS	Description	
Warranty	Limited Lifetime	