

## Non-Tilt Pneumatic Cylinder Overbed Table U

HCPCS:



- Pneumatic cylinder lift mechanism infinite height adjustment from 27"42"
- ¾" x 18" x 32" high pressure laminate top with T- molding
- 2" casters for easy maneuvering
- Column can be mounted on left or right side
- Optional dual mirror vanity drawer
- Available with "H" (Figure A) or "U" base in Chrome or Opal powder coat
- Base and columns available in additional powder coat finishes
- Available on drop ship basis only

| PRODUCTS    | Description                     |
|-------------|---------------------------------|
| 16008-UCBO  | OBT-PnmCyl-BO-U-Ch-SpclOrdr     |
| 16008-UCBOV | OBT-PnmCyl-BO-U-Ch-Van-SpclOrdr |
| 16008-UCDM  | OBT-PnmCyl-FM-U-Ch-SpclOrdr     |
| 16008-UCDMV | OBT-PnmCyl-FM-U-Ch-Van-SpclOrdr |
| 16008-UCEO  | OBT-PnmCyl-EO-U-Ch-SpclOrdr     |
| 16008-UCEOV | OBT-PnmCyl-EO-U-Ch-Van-SpclOrdr |
| 16008-UCFM  | OBT-PnmCyl-FM-U-Ch-SpclOrdr     |
| 16008-UCFMV | OBT-PnmCyl-FM-U-Ch-Van-SpclOrdr |
| 16008-UCGW  | OBT-PnmCyl-GW-U-Ch-SpclOrdr     |
| 16008-UCGWV | OBT-PnmCyl-GW-U-Ch-Van-SpclOrdr |
| 16008-UCNT  | OBT-PnmCyl-NT-U-Ch-SpclOrdr     |
| 16008-UCNTV | OBT-PnmCyl-NT-U-Ch-Van-SpclOrdr |
| 16008-UCQO  | OBT-PnmCyl-QO-U-Ch-SpclOrdr     |
| 16008-UCQOV | OBT-PnmCyl-QO-U-Ch-Van-SpclOrdr |
| 16008-UCSO  | OBT-PnmCyl-SO-U-Ch-SpclOrdr     |
| 16008-UCSOV | OBT-PnmCyl-SO-U-Ch-Van-SpclOrdr |
| 16008-UCWC  | OBT-PnmCyl-WC-U-Ch-SpclOrdr     |
| 16008-UCWCV | OBT-PnmCyl-WC-U-Ch-Van-SpclOrdr |
| 16008-UOBO  | OBT-PnmCyl-BO-U-Op-SpclOrdr     |
| 16008-UOBOV | OBT-PnmCyl-BO-U-Op-Van-SpclOrdr |
| 16008-UODM  | OBT-PnmCyl-FM-U-Op-SpclOrdr     |
| 16008-UODMV | OBT-PnmCyl-FM-U-Op-Van-SpclOrdr |
| 16008-UOEO  | OBT-PnmCyl-EO-U-Op-SpclOrdr     |
| 16008-UOEOV | OBT-PnmCyl-EO-U-Op-Van-SpclOrdr |
| 16008-UOFM  | OBT-PnmCyl-FM-U-Op-SpclOrd      |

## PRODUCTS

## Description

|             |                                  |
|-------------|----------------------------------|
| 16008-U0FMV | OBT-PnmCyl-FM-U-Op-Van-SpclOrd   |
| 16008-U0GW  | OBT-PnmCyl-GW-U-Opal             |
| 16008-U0GWV | OBT-PnmCyl-GW-U-Opal-Van         |
| 16008-U0NT  | OBT-PnmCyl-NT-U-Op-SpclOrdr      |
| 16008-U0NTV | OBT-PnmCyl-NT-U-Op-Van-SpclOrdr  |
| 16008-U0QO  | OBT-PnmCyl-QO-U-Op-SpclOrdr      |
| 16008-U0QOV | OBT-PnmCyl-QO-U-Op-Van-SpclOrdr  |
| 16008-U0SO  | OBT-PnmCyl-SO-U-Op-SpclOrder     |
| 16008-U0SOV | OBT-PnmCyl-SO-U-Op-Van-SpclOrder |
| 16008-U0WC  | OBT-PnmCyl-WC-U-Op-SpclOrdr      |
| 16008-U0WCV | OBT-PnmCyl-WC-U-Op-Van-SpclOrdr  |