

## Table de lit à vérin pneumatique non inclinable

### U

HCPCS:



- Mécanisme de levage à vérin pneumatique réglage de hauteur infini de 27" à 42"
- Dessus en stratifié haute pression ¾" x 18" x 32" avec moulure en T
- Roulettes de 2" pour une manuvre facile
- La colonne peut être montée sur le côté gauche ou droit
- Disponible avec une base en "H" (Figure A) ou en "U" en chrome ou en poudre opale < /li>
- Tiroir de meuble-lavabo à double miroir en option
- Base et colonnes disponibles dans des finitions de revêtement en poudre supplémentaires
- Disponible uniquement en livraison directe

| PRODUCTS    | Description                     |
|-------------|---------------------------------|
| 16008-UCBO  | OBT-PnmCyl-BO-U-Ch-SpclOrdr     |
| 16008-UCBOV | OBT-PnmCyl-BO-U-Ch-Van-SpclOrdr |
| 16008-UCDM  | OBT-PnmCyl-FM-U-Ch-SpclOrdr     |
| 16008-UCDMV | OBT-PnmCyl-FM-U-Ch-Van-SpclOrdr |
| 16008-UCEO  | OBT-PnmCyl-EO-U-Ch-SpclOrdr     |
| 16008-UCEOV | OBT-PnmCyl-EO-U-Ch-Van-SpclOrdr |
| 16008-UCFM  | OBT-PnmCyl-FM-U-Ch-SpclOrdr     |
| 16008-UCFMV | OBT-PnmCyl-FM-U-Ch-Van-SpclOrdr |
| 16008-UCGW  | OBT-PnmCyl-GW-U-Ch-SpclOrdr     |
| 16008-UCGWV | OBT-PnmCyl-GW-U-Ch-Van-SpclOrdr |
| 16008-UCNT  | OBT-PnmCyl-NT-U-Ch-SpclOrdr     |
| 16008-UCNTV | OBT-PnmCyl-NT-U-Ch-Van-SpclOrdr |
| 16008-UCQO  | OBT-PnmCyl-QO-U-Ch-SpclOrdr     |
| 16008-UCQOV | OBT-PnmCyl-QO-U-Ch-Van-SpclOrdr |
| 16008-UCSO  | OBT-PnmCyl-SO-U-Ch-SpclOrdr     |
| 16008-UCSOV | OBT-PnmCyl-SO-U-Ch-Van-SpclOrdr |
| 16008-UCWC  | OBT-PnmCyl-WC-U-Ch-SpclOrdr     |
| 16008-UCWCV | OBT-PnmCyl-WC-U-Ch-Van-SpclOrdr |
| 16008-UOBO  | OBT-PnmCyl-BO-U-Op-SpclOrdr     |
| 16008-UOBOV | OBT-PnmCyl-BO-U-Op-Van-SpclOrdr |
| 16008-UODM  | OBT-PnmCyl-FM-U-Op-SpclOrdr     |
| 16008-UODMV | OBT-PnmCyl-FM-U-Op-Van-SpclOrdr |
| 16008-UOEO  | OBT-PnmCyl-EO-U-Op-SpclOrdr     |
| 16008-UOEOV | OBT-PnmCyl-EO-U-Op-Van-SpclOrdr |
| 16008-UOFM  | OBT-PnmCyl-FM-U-Op-SpclOrd      |

## PRODUCTS

## Description

|             |                                      |
|-------------|--------------------------------------|
| 16008-U0FMV | OBT-PnmCyl-FM-U-Op-Van-SpclOrd       |
| 16008-U0GW  | OBT-PnmCyl-GW-U-Opale                |
| 16008-U0GWV | OBT-PnmCyl-GW-U-Opal-Van             |
| 16008-U0NT  | OBT-PnmCyl-NT-U-Op-SpclOrd           |
| 16008-U0NTV | OBT-PnmCyl-NT-U-Op-Van-SpclOrd       |
| 16008-U0QO  | OBT-PnmCyl-QO-U-Op-SpclOrd           |
| 16008-U0QOV | OBT-PnmCyl-QO-U-Op-Van-SpclOrd       |
| 16008-U0SO  | Commande OBT-PnmCyl-SO-U-Op-Spcl     |
| 16008-U0SOV | Commande OBT-PnmCyl-SO-U-Op-Van-Spcl |
| 16008-U0WC  | OBT-PnmCyl-WC-U-Op-SpclOrd           |
| 16008-U0WCV | OBT-PnmCyl-WC-U-Op-Van-SpclOrd       |