



COLORADO SPRING BATHING SYSTEM

The Colorado Springs Bathing System is designed to provide a unique combination of function, comfort, and luxury. It is a freestanding tub that allows the caregiver 360° unobstructed access to the resident. The low threshold and wide side-entry door allows easy access for residents with various degrees of mobility and the built-in contoured seat allows the resident to sit comfortably and safely in a slightly reclined position. The smooth and soft lines will enhance the look of the bathing area and allow for easy maintenance and cleaning. Designed to create the atmosphere of a luxury spa, the sleek exterior design of the Colorado Springs Bathing System will give the bathing area an upscale look. Caregivers will enjoy the Colorado Springs Bathing System for its ease of use and the improved quality of care they can provide their residents.

**600 lb.
WEIGHT CAPACITY**

DB-COLORADO-1000	Side-Entry with Thermostatic Mixing Valve
DB-COLORADO-1100	Side-Entry with Thermostatic Mixing Valve and Disinfectant System
DB-COLORADO-1200	Side-Entry with Thermostatic Mixing Valve, Disinfectant System and Air Spa

FEATURES:

- ◆ **Thermostatic mixing valve** - Provides an anti-scald protection that dispenses 15 gallons per minute for fast tub fill.
- ◆ **Built-in disinfection system** - The most time-efficient solution for disinfection of the bath's external surfaces and minimizes the risk of cross-infection. The built-in disinfection system cabinet is conveniently located on the front of the tub for easy access. The lockable cabinet provides a secure place for cleaning supplies and disinfectants.
- ◆ **Shower wand** - Hair and upper body washing can be used at the same time the tub is filling.
- ◆ **Water temperature display** - Monitors tub-fill and shower wand temperatures.
- ◆ **3/4" supply lines** - Constructed of stainless steel mesh clad.
- ◆ **Lap and shoulder harness** - Used for resident safety.
- ◆ **Secondary door lock** - Prevents accidental door opening when in use.
- ◆ **The CleanDrive Air Spa System** - Variable speeds and pulse modes designed to enhance the bathing experience for the resident and caregiver.

Additional features of the CleanDrive Air System include:

- ◆ The hydro massage provides millions of preheated air bubbles in the water which increases stimulation of peripheral blood circulation, reduces pain and inflammation in the joints, gently massages away dirt, and induces deeper relaxation.
- ◆ The CleanDrive Air Spa System is a hygienic and low maintenance hydro massage system featuring 12 air jets with spring loaded check valves for back flow prevention. It prevents water from entering the air lines, eliminating the extra step to disinfect air lines. The CleanDrive Air Spa System simplifies disinfection and reduces the possibility of cross-contamination.
- ◆ Unlike most other whirlpool and air spa systems, the CleanDrive Air Spa System does not allow water to enter the air lines to reduce the growth of bacteria and mold.
- ◆ The motor is quiet, efficient, and has a ceramic heating element that blows warm air to help maintain the water temperature and warm the surface of the tub.
- ◆ Warranty: 5 Years on finishes and composite shell, 3 years on all moving and stationary parts, and lifetime on the door seal.

To order, call LTCCare Customer Service Toll Free at (877) 374-8763, contact your Sales Representative or visit us at www.drivemedical.com TODAY!



The built-in disinfection system cabinet is conveniently located on the side of the tub for easy access. The lockable cabinet provides a secure place for cleaning supplies and disinfectants.



The wide side-entry door and unobstructed mobile lift access will allow the resident to enter and exit the bathtub easily and safely.



The Colorado Spring features an easy-to-use control panel which includes: the tub-fill and shower valves, water temperature display, thermostatic mixing valve, and CleanDrive Air Spa System control pad.



The contoured seat, lap and shoulder harness, and strategically placed air spa jets will provide a safe and pleasurable experience for your resident.

**Made
in the
USA**

SPECIFICATIONS

Width:	33.5"
Length:	61.5"
Height (to top of console):	42.5"
Height (to top of door):	35"
Seat Height:	21" (from floor)
Seat Width:	21.5"
Threshold Height:	11"
Transfer Device Access:	7" H x 44" W
Water Supply Inlet:	3/4" diameter
Drain Size:	1.5"
Water Capacity (Occupied):	48 gallons
Fill Time:	3 minutes, 12 seconds
Drain Time:	1 minute, 58 seconds
Weight Capacity (resident):	600 lbs
Door Configuration:	Right Hinge only
Electrical Supply:	115V AC; 60Hz; 15 amp

To order, call LTCCare Customer Service Toll Free at (877) 374-8763, contact your Sales Representative or visit us at www.drivemedical.com TODAY!