



SARATOGA SPRING BATHING SYSTEM

The Saratoga Springs Bathing System is designed for care facilities looking for a more time efficient and labor-saving bathing solution. It has a unique tilting function, which allows patients with challenging conditions to recline in a contoured seat and become immersed in warm water. At only 24 gallons to pre-fill the foot well, water usage and the associated costs are greatly reduced. The sleek exterior design of the Saratoga Springs Bathing System will give the bathing area the look of a luxury spa and allows for easy maintenance and cleaning. Caregivers will enjoy the Saratoga Springs Bathing System for its ease of use and the improved quality of care they can provide their larger residents.

**450 lb.
WEIGHT CAPACITY**

DB-SARATOGA-1000	Reclining Tub with Thermostatic Mixing Valve and Auto-Fill
DB-SARATOGA-1100	Reclining Tub with Thermostatic Mixing Valve, Disinfectant System and Auto-Fill
DB-SARATOGA-1200	Reclining Tub with Thermostatic Mixing Valve, Disinfectant System, Air Spa and Auto-Fill

FEATURES:

- ◆ **Auto-fill** - Eliminates waiting and/or observation time when pre-filling the bathtub to a pre-set level.
- ◆ **Thermostatic mixing valve** - Provides an anti-scald protection that dispenses 15 gallons per minute for fast tub fill.
- ◆ **Built-in disinfection system** - The most time-efficient solution for disinfection of the bath's external surfaces and minimizes the risk of cross-infection. The built-in disinfection system cabinet is conveniently located on the back end of the tub for easy access. The lockable cabinet provides a secure place for cleaning supplies and disinfectants.
- ◆ **Handheld shower wand** - Hair and upper body washing can be used at the same time the tub is filling.
- ◆ **Water temperature display** - Monitors tub-fill and shower wand temperatures.
- ◆ **3/4" supply lines** - Constructed of stainless steel mesh clad.
- ◆ **Three plumbing configurations** - to choose from side, front, or island configuration. The island configuration allows 360° unobstructed access for the caregiver and with no visible hook-ups or cords, the Saratoga Springs Bathing System will give the bathing area the look of a luxury spa.
- ◆ **The CleanDrive Air Spa System** - Variable speeds and pulse modes designed to enhance the bathing experience for the resident and caregiver.

Additional features of the CleanDrive Air System include:

- ◆ The hydro massage provides millions of preheated air bubbles in the water which increases stimulation of peripheral blood circulation, reduces pain and inflammation in the joints, gently massages away dirt, and induces deeper relaxation.
- ◆ The CleanDrive Air Spa System is a hygienic and low maintenance hydro massage system featuring 12 air jets with spring loaded check valves for back flow prevention. It prevents water from entering the air lines, eliminating the extra step to disinfect air lines. The CleanDrive Air Spa System simplifies disinfection and reduces the possibility of cross-contamination.
- ◆ Unlike most other whirlpool and air spa systems, the CleanDrive Air Spa System does not allow water to enter the air lines to reduce the growth of bacteria and mold.
- ◆ The motor is quiet, efficient, and has a ceramic heating element that blows warm air to help maintain the water temperature and warm the surface of the tub.
- ◆ Warranty: 5 Years on finishes and composite shell, 3 years on all moving and stationary parts, and lifetime on the door seal.

To order, call LTCCare Customer Service Toll Free at (877) 374-8763, contact your Sales Representative or visit us at www.drivemedical.com TODAY!



The Saratoga Spring features an easy-to-use control panel which includes: the tub-fill and shower valves, water temperature display, thermostatic mixing valve, and CleanDrive Air Spa System control pad.



The Saratoga Spring has a unique tilting function, which allows patients with challenging conditions to recline in a contoured seat and become immersed in warm water.



The sleek design of the Saratoga Spring will give your bathing area the look of a luxury spa.



The built-in disinfection system cabinet is conveniently located on the end of the tub for easy access. The lockable cabinet provides a secure place for cleaning supplies and disinfectants.

Made
in the
USA

SPECIFICATIONS

Width:	34"
Length:	73"
Length Reclined:	88"
Height:	51"
Height with Door Open:	78"
Seat Height:	21" (from floor)
Seat Width:	24"
Transfer Device Access:	7" H x 31" W
Water Supply Inlet:	3/4" diameter
Drain Size:	1.5"
Water Capacity (occupied):	38 gallons
Water Capacity (foot well):	24 gallons
Fill Time (foot well):	1 minute, 36 seconds
Drain Time (foot well):	2 minutes, 10 seconds
Weight Capacity (resident):	450 lbs
Weight Capacity (resident & water):	730 lbs
Door Configuration:	Upward Swing, Left Hinge only
Electrical Supply:	115V AC; 60Hz; 15 amp

To order, call LTCCare Customer Service Toll Free at (877) 374-8763, contact your Sales Representative or visit us at www.drivemedical.com TODAY!