


SCOUT 3



1. Speed Dial- BC31001
Speed Dial Knob- LRM412103-11
2. Throttle Assembly Old- [2]
Throttle Assembly New- S12-069-00300
3. Battery Gauge (4 LEDs)- [5]
Battery Gauge (8 LEDs)- C11-076-00100 [1]
4. Horn Button- SC31034
5. Key Switch Assembly- BC31006
Key- BCKEY
6. Upper Tiller Cover- SC31008
7. Lower Tiller Cover- SC31009
8. Basket Bracket- BC35051
9. Basket- S35003-P
10. Handgrip Foam- SC31002
Handgrip Plastic- SC31054
11. Tiller- SC31001
12. Tiller Swivel Assembly- SC31010
13. Tiller Lock- C08-069-01000
14. Tiller Bearing Set- BC31021
15. Fork- SC31038
16. Front Wheel Black- SC31043
Deluxe Front Wheel- SC31044-D
17. Charger- LRM302111
18. Front Axle- SC31046
19. Front Shroud- SC31039
20. Foot Shroud- SC31024
21. Battery Box Shroud- SC31027
22. Disconnect Handle Assembly- SC31019
23. Red Decoration- SC31041
Deluxe Red Decoration- SC31041-D
Blue Decoration- SC31042
Deluxe Blue Decoration- SC31042-D
24. Front Frame- SC31040
25. Black Seat- SC35045
Deluxe Two-Tone Seat- SC31035
Deluxe Back Cover- [5]
Deluxe Seat Cover- [5]
26. Left Armrest- BC35040L
Right Armrest- BC35040R
Arm Pad- BC953014
27. Armrest Knob- S31052
28. Seat Plate- SC31023
29. Seat Post- SC31018 [4]
30. Seat Post Bolt- C07-071-00101
31. Battery Box Hold Down- SC31047
32. Battery Box & Batteries (12V 12ah)- SC31007
Battery Box Only (12V 12ah)- SC31006
Battery 12V 12ah (2 needed)- LRM402109
EXT model
Battery Box & Batteries (12V 20ah)- SC31050
Battery Box Only (12V 20ah)- SC31051
Battery 12V 20ah (2 needed)- LRM412123
33. Charger Socket & Harness- SC31031 [6]
34. 20 Amp Circuit Breaker- BC20ACB
35. Battery Contact- C12-060-00300
36. Motor/Brake/Transaxle Assembly- SC31053
Transaxle- BC350112
Motor- SC31036
Electro Magnetic Brake- SC31037
37. Rear Wheel Black- SC31003
Deluxe Rear Wheel- SC31003-D
38. Hub Cap- SC31025
39. Rear Axle Bolt Assembly- SC31049
40. Controller Curtis- SC31052
Controller Dynamic- C10-069-00100-3 [3]
41. Wire Harness Complete-
Curtis- SC31026
Dynamic- DST-018 [3]
42. Contact Block Set w/Motor Harness- SC31033
43. Rear Shroud- SC31011
44. Rear Frame- SC31012
45. Anti-Tipper Wheel- S31026
46. Motor Wire Harness- SC31020

[1] After Serial # VP035868

[2] Stock is gone use new style.

[3] After Serial # VP142824

[4] No longer available in chrome.

[5] Part no longer available.

[6] For 12ah battery box only.