

Cobra GT4 4-Wheel Heavy Duty Scooter

Order Form/Price List



Company Name:		Phone:
Contact Person:	P.O. #	Date:
Bill to:		Ship to:

Scooter with Seating

Captain's Seat

- ☐ COBRAGT420CS..... 20" (W) x 20" (D)\$8,645.00
- ☐ COBRAGT422CS..... 22" (W) x 20" (D)\$9,009.00

Batteries & Charger

- ☐ 12V x 75AH Batteries (pair)Included
- ☐ 5A Battery ChargerIncluded

Accessories

- ☐ AB1000 Armrest Bag\$25.00
- ☐ AB1120 Backpack\$56.00
- ☐ AH1000 Crutch/Cane Holder\$53.00
- ☐ AH2000 Cylinder Tank Carrier\$70.00
- ☐ AB2100 Rear Basket\$70.00
- ☐ AZ4000 Scooter Cover\$72.00
- ☐ ST301-COVER 20" Seat Cover\$118.17
- ☐ ST306-COVER 22" Seat Cover\$118.17
- ☐ AT1000 Scooter Trailer\$204.00
- ☐ AZ0060 Cup Holder\$30.00

Specifications

HCPCS	N/A
Performance	
Maximum Weight Capacity	450 lbs.
Maximum Speed	10 mph
Maximum Range	35 miles
Maximum Grade/Incline	7.5°
Turning Radius	82"
Dimensions	
Overall Length	63"
Overall Width	25"
Floor to Seat Height	27"-30"
Ground Clearance	7"
Power	
Motor	24V x 950W x 5200rpm
Controller	Dynamic Rhino II 160A
Brakes	Electromagnetic
Freewheel Mode	Yes
Tires (Pneumatic)	
Front	14" x 4"
Rear	14" x 4"
Weights	
Base Weight	191 lbs.
Battery Weight	94 lbs.
Seat Weight	55 lbs.
Total Weight	340 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.