


770-525 AQUASENSE FOLDING BATH SEAT WITH BACK



Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	B04-252	Push Button & Stainless Steel Spring (4 / set)	1	1	ST
2.	A04-265	Anti-Rattle, Blue (4 / set)	1	1	ST
3.	A04-401	Rubber Tip, Blue (4 / set)	1	1	ST
4.	B04-001	Extension Leg w/Hardware (Aluminum Tube ϕ 28.7x255mm, Blue Anti-Rattle & Blue Tip) (4 / set)	1	1	ST
5.	B04-101	Seat Retaining Clip w/Hardware (1 Clip & 2 Self-Threading Screws) (4 / set)	1	1	ST