



PreserveTech™
Stays Cleaner • Lasts Longer

Raised Toilet Seat with Bidet

The PreserveTech™ Raised Toilet Seat with Bidet offers 4" of added height to make transfers on and off the toilet more accessible and comfortable. The bidet provides front and back cleansing, and the control panel allows the user to adjust water pressure, cleansing modes, nozzle cleaning, and water temperature (RTL12C006-WH only). The contoured seat fits most standard and elongated toilets and includes a lid for added hygiene and aesthetic.

EASY-TO-USE CONTROLS
Activate and adjust water pressure and select cleansing modes

HYGIENIC & GENTLE
Provides front or back cleansing depending on needs

EXTRA HEIGHT
Contoured seat adds 4" of height to standard or elongated toilets for safer transfer

FITS MOST STANDARD & ELONGATED TOILETS

Drive PreserveTech™ products are manufactured with an antimicrobial* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.

*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

USER-FRIENDLY
Easily connects to water source*, no electric connection or batteries required



| ITEM # | PRODUCT DESCRIPTION | UPC | UOM |
|--------------|---|-------------------------------------|------------|
| RTL12C005-WH | PreserveTech™ Raised Toilet Seat with Bidet, Ambient Water Temperature | 822383015941/ea / 50822383015946/cv | 1/ea, 1/cv |
| RTL12C006-WH | PreserveTech™ Raised Toilet Seat with Bidet, Ambient & Warm Water Temperature | 822383015958/ea / 50822383015953/cv | 1/ea, 1/cv |

PRODUCT SPECS

| | |
|--------------------|--|
| Product Dimensions | 17.14" (W) x 16.13" (L) x 7.58" (H) |
| Product Weight | 3.5 lb. (RTL12C005-WH) 3.54 lb. (RTL12C006-WH) |
| Color | White |
| Weight Capacity | 400 lb. |
| Warranty | One Year Limited** |

BOX SPECS

| | |
|-----------------------|--|
| Retail Box Dimensions | 8.07" (W) x 17.52" (L) x 16.54" (H) |
| Retail Box Weight | 4.5 lb. (RTL12C005-WH) 4.54 lb. (RTL12C006-WH) |

*Initial temperature of water release is contingent upon current temperature of surroundings and will warm based upon your current water setting and supply.

**Refer to instruction guide or website for warranty details.



SCAN FOR BILT. 3D INSTRUCTIONS



POWERED BY BILT INTELLIGENT INSTRUCTIONS*

FOR USE WITH ITEM #RTL12C005-WH

SCAN FOR BILT. 3D INSTRUCTIONS



POWERED BY BILT INTELLIGENT INSTRUCTIONS*

FOR USE WITH ITEM #RTL12C006-WH

