



## 770-500 AQUASENSE ADJUSTABLE BATH BENCH




Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	<b>A11-401</b>	Seat, Blow-Molded Plastic, White (505x305x35mm)	1	1	EA
2.	<b>B11-822</b>	Hardware Set, Aqua-Blue (5 Triangular Knobs M6x32mm)	1	1	ST
3.	<b>B11-811</b>	Push Button dia:8mm & Spring - 4ea / set	1	1	ST
4.	<b>A11-212</b>	Anti-Rattle for dia:28.6mm, Aqua-Blue – 4 ea / set	1	1	ST
5.	<b>B11-012</b>	Extension Leg, Aqua-Blue (Aluminum Tube dia:28.6x250mm, Anti-Rattle & Rubber Tip) – 4 ea / set	1	1	ST
6.	<b>A11-482</b>	Rubber Tip for dia:28.6mm, Aqua-Blue – 4 ea / set	1	1	ST
Not Shown	<b>B11-851</b>	Hardware Set for Bath Bench (5 Cross-Cap Screws M6x32mm, 5 Curved Washers & 1 Locknut) -- <b>Accessory (screw set to replace triangular knobs)</b>	1	1	ST

## 770-501 ADJUSTABLE BATH BENCH




Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	<b>A03-452</b>	Seat / Backrest, Blow-Molded Plastic, White (505x305x35mm)	1	1	EA
2.	<b>B03-815</b>	Hardware for Legs (Bolt M6x37mm, Tooth-Lock Washer, S.Steel Curved Washer & Nut)	1	1	ST
3.	<b>B03-011</b>	Top Leg (Aluminum Tube $\phi$ 25.4mm & 2 Push-Buttons with 2 Springs)	1	1	EA
4.	<b>B03-012</b>	Bottom Leg (Aluminum Tube $\phi$ 25.4mm & 2 Push-Buttons with 2 Springs)	1	1	EA
5.	<b>B03-816</b>	Hardware for 770-500 Seat / 770-510 Back (Screw 1/4"x38mm & S.Steel Curved Washer)	4	1	ST
6.	<b>B03-811</b>	Push Button & Spring	4	1	ST
7.	<b>A03-210</b>	Anti-Rattle for Extension Leg, Grey Plastic	4	1	EA
8.	<b>A03-432</b>	Rubber Tip, Grey	4	1	EA
9.	<b>B03-010</b>	Leg Extension Assy (Aluminum Tube $\phi$ 28.6, Anti-Rattle & Rubber Tip)	4	1	EA
Not shown	<b>B03-821</b>	Hardware Set for 770-500 (Bolt M6x37, 4 Screws 1/4"x38mm, Tooth Lock Washer, 5 S.Steel Curved Washers & Nut)	1	1	ST

## 770-506 AQUASENSE ADJUSTABLE BATH BENCH - BLUE




Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	<b>A11-402</b>	Seat, Blow-Molded Plastic, Aqua-Blue (505x305x35mm)	1	1	EA
2.	<b>B11-822</b>	Hardware Set, Aqua-Blue (5 Triangular Knobs M6x34mm)	1	1	ST
3.	<b>B11-811</b>	Push Button dia:8mm & Spring - 4ea / set	1	1	ST
4.	<b>A11-212</b>	Anti-Rattle for dia:28.6mm, Aqua-Blue - 4 ea / set	1	1	ST
5.	<b>B11-012</b>	Extension Leg, Aqua-Blue (Aluminum Tube dia:28.6x250mm, Anti-Rattle & Rubber Tip) - 4 ea / set	1	1	ST
6.	<b>A11-482</b>	Rubber Tip for dia:28.6mm, Aqua-Blue - 4 ea / set	1	1	ST
Not Shown	<b>B11-851</b>	Hardware Set for Bath Bench w/out Back (5 Cross-Cup Screws M6x34mm, 5 Curved Washers) - <b>Accessory</b>	1	1	ST

## 770-507 AQUASENSE ADJUSTABLE BATH BENCH - YELLOW




Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	<b>A11-403</b>	Seat, Blow-Molded Plastic, Yellow (505x305x35mm)	1	1	EA
2.	<b>B11-823</b>	Hardware Set, Yellow (5 Triangular Knobs M6x34mm)	1	1	ST
3.	<b>B11-811</b>	Push Button dia:8mm & Spring - 4ea / set	1	1	ST
4.	<b>A11-213</b>	Anti-Rattle for dia:28.6mm, Yellow - 4 ea / set	1	1	ST
5.	<b>B11-013</b>	Extension Leg, Yellow (Aluminum Tube dia:28.6x250mm, Anti-Rattle & Rubber Tip) - 4 ea / set	1	1	ST
6.	<b>A11-483</b>	Rubber Tip for dia:28.6mm, Yellow - 4 ea / set	1	1	ST
Not Shown	<b>B11-851</b>	Hardware Set for Bath Bench w/out Back (5 Cross-Cup Screws M6x34mm, 5 Curved Washers) - <b>Accessory</b>	1	1	ST

## 770-508 AQUASENSE ADJUSTABLE BATH BENCH - TAUPE



Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	<b>A11-404</b>	Seat, Blow-Molded Plastic, Taupe (505x305x35mm)	1	1	EA
2.	<b>B11-824</b>	Hardware Set, Taupe (5 Triangular Knobs M6x34mm)	1	1	ST
3.	<b>B11-811</b>	Push Button dia:8mm & Spring - 4ea / set	1	1	ST
4.	<b>A11-214</b>	Anti-Rattle for dia:28.6mm, Taupe - 4 ea / set	1	1	ST
5.	<b>B11-014</b>	Extension Leg, Taupe (Aluminum Tube dia:28.6x250mm, Anti-Rattle & Rubber Tip) - 4 ea / set	1	1	ST
6.	<b>A11-484</b>	Rubber Tip for dia:28.6mm, Taupe - 4 ea / set	1	1	ST
Not Shown	<b>B11-851</b>	Hardware Set for Bath Bench w/out Back (5 Cross-Cup Screws M6x34mm, 5 Curved Washers) - <b>Accessory</b>	1	1	ST

## 770-509 AQUASENSE ADJUSTABLE BATH BENCH



Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	<b>A11-401</b>	Seat, Blow-Molded Plastic, White (505x305x35mm)	1	1	EA
2.	<b>B11-822</b>	Hardware Set, Aqua-Blue (5 Triangular Knobs M6x32mm)	1	1	ST
3.	<b>B11-811</b>	Push Button dia:8mm & Spring - 4ea / set	1	1	ST
4.	<b>A11-212</b>	Anti-Rattle for dia:28.6mm, Aqua-Blue – 4 ea / set	1	1	ST
5.	<b>B11-012</b>	Extension Leg, Aqua-Blue (Aluminum Tube dia:28.6x250mm, Anti-Rattle & Rubber Tip) – 4 ea / set	1	1	ST
6.	<b>A11-482</b>	Rubber Tip for dia:28.6mm, Aqua-Blue – 4 ea / set	1	1	ST
Not Shown	<b>B11-851</b>	Hardware Set for Bath Bench (5 Cross-Cap Screws M6x32mm, 5 Curved Washers & 1 Locknut) -- <b>Accessory (screw set to replace triangular knobs)</b>	1	1	ST