


STSP450



- | | |
|---|---|
| 1. Upper Mast Assy.- STS450-N-10001AE | 12. Control Box- STS450-N-15E |
| 2. Foot Plate- STS450-N-10004AE | 13. Hand Control- STS450-N-20E |
| 3. Knee Pad Bracket Assy.- STS450-N-10005AE | 14. Charger Stand- STS450-N-21E |
| 4. Base Assy. Left- STS450-N-10006E | 15. External Charger- STS450-N-22E |
| 5. Base Assy. Right- STS450-N-10007E | 16. Power Base Cord- STS450-N-23E |
| 6. Knee Pad Assy.- STS450-N-10008E | 17. Leg Strap Bracket- STS450-N-10024E |
| 7. Linkage Connection plate- STS450-N-10009E | 18. Non-Locking Caste- UA17580750E |
| 8. Foot Actuator Assy.- STS450-N-10016E | 19. Locking Caster- UA18661000E |
| 9. Elevation Actuator Assy.- STS450-N-10017AE | 20. Battery- SP07-SP-2139-002 |
| 10. Base Actuating Rod- STS450-N-11E | 21. Receiver Plate Assy.- STS450-N-10003E |
| 11. Connection Linkage Bushing- STS450-N-13E | Not Shown Leg Strap- STSSTRAP |