


301



		2M or 20 H Serial #	3H Serial #
1.	Hand Brake	9502M30103	10208BRAKE
2.	Hand Grips (pr)	9501S72602	[1]
3.	Height Adj. Knob	9502M30105	10208KN
4.	Rear Wheel	9502M30114	[1]
5.	Seat & Brackets	9502M30115	[1]
6.	Seat Brackets	9505W1026104	9505W1026104
7.	a. Spring Boot	[1]	[1]
	b. Adjuster	9500001232	9500001232
	c. Fastener	9500004101	9500004101
8.	Front Wheel	9502M30113	[1]
9.	Front Wheel Bearing	9502M30111	STDS1011
10.	Fork	9502M30112	[1]
11.	Fork Stem Bearing	STDS1060	STDS1060
12.	Backrest	9502M30116	[1]
13.	Anti-Rattler	950300169901	950300169901
14.	Cable & Sleeve	[1]	950300089901
15.	Front Axle	9502M30110	STDS2A3423
16.	Fork Fastener	9502M30119	[1]
17.	Rear Axle	9502M30109	[1]

[1] Stock is gone the part is no longer available.