


FLP600



- 1. 3" Caster- SP07-2AAX0134
- 2. 5" Locking Caster- SP07-2AAX0133
- 3. Boom Motor- SP07-TA12-2139-001
- 4. Leg Spreader Motor- SP07-TA16-2139-001
- 5. Controller- SP07-TC12-2139-002
- 6. Hand Pendant- SP07-TH10-2139-001
- 7. Battery 5.0 Ah- SP07-SP-2139-002
- 8. Charger- SP07-TBC2-42-211-1
- 9. Charger Power Cord- SP07-TBC2-2139-001
- 10. Mast/Boom Assembly Cover- SP07-2MBA0009
- 11. Base Cover Bolt Set (4 ea.)- SP07-2BCB0008

- 12. Push Handle with Hardware- SP07-1L30080F
- 13. Cradle Pivot- 13058PB-FL
- 13. Cradle Pivot Attachment Bolt- SP07-350F0190-49
- 14. Cradle Holder Hardware- SP07-350F0554
- 15. Motor to Controller Y Cable- SP07-2LTE0238
- 16. Base w/o Wheels- SP07-FLPBASE [1]
- 17. Boom- SP07PUB
- 18. Mast- SP07PLM
- 19. Actuator Bolt Set (2 ea)- SP07-2139-001BOLTS
- 20. Cradle- 13057
- 21. Controller Mount- 13260-TM

[1] Includes motor, part 4.