

Disposable Container with Internal Filter



(Exploded View)



Assembled Collection Container Kit 7305D-633 - Includes internal filter cartridge, 800 ml container and lid, splash guard, and 4.5" connection and 6' patient tubing (1 ea)

Assembled Collection Container Kit 7305D-633-25 - Includes internal filter cartridge, 800 ml container and lid, splash guard, and 4.5" connection and 6' patient tubing (25 pk)



Unassembled Collection Container 7305D-632 - Includes internal filter cartridge, 800 ml container and lid, splash guard, and 4.5" connection tubing (48 pk)



Elbow with Integrated Filter 7305D-635 (12/pk)



6' Patient Tubing SUCP TUBING 72 (1/ea)

Disposable Container with External Filter



(Exploded View)



Assembled Collection Container Kit 22330 - Includes 800 ml disposable container & lid, external bacteria filter, elbow, 10" connection and 6' patient tubing package (1/ea)

Assembled Collection Container Kit 22330-12 - 800 ml disposable container & lid, external bacteria filter, elbow, 10" connection and 6' patient tubing package (12/pk)



Disposable Container and Lid only 610-12B (12/pk)
610-48BP (48/pk)



6' Patient Tubing SUCP TUBING 72 (1/ea)



Universal Suction Tubing and Filter Kit 18600-KITN - (filter, elbow, 10" connection and 6' patient tubing package) (1/ea)



3.375" Connection Tubing 7305D-612 (6/pk)



Filter 18600-FILTER (3/pk)



Filter 7305D-608 (12/pk)



Elbow 7305D-619 (6/pk)



Elbow, Filter, 3.375" Connection Tubing Kit 7305D-614 (1/ea)

NOTE - The 18600 Drive Heavy-Duty Suction Unit requires the 10" tubing connection included in the 22330 container series only. Any container series may be used with the 7325, 7305 and 7314 DeVilbiss Vacu-Aide® Suction Units.

Accessories and Replacement Parts



7325 Carrying Case
7325D-635



7314 Carrying Case
7314D-606



7305 Carrying Case
7305D-606



12V DC Power Cord
7304D-619



AC to DC Power Adapter/Charger
7314P-613



Power Cord for USA
DV51D-606



Vacuum Gauge Test Fixture
7304D-636



Connection Tubing
SUCP TUBING 10

Battery with Foam for 7314P-D
7314P-614

Battery for 7305P-D
7305P-614

Battery for 7325P-D
7325P-614

Battery Retrofit Kit for 7305D-D
7305P-620