



**PreserveTech™**  
Stays Cleaner • Lasts Longer

# Raised Toilet Seat with Bidet

The PreserveTech™ Raised Toilet Seat with Bidet offers 4" of added height to make transfers on and off the toilet more accessible and comfortable. The bidet provides front and back cleansing, and the control panel allows the user to adjust water pressure, cleansing modes, nozzle cleaning, and water temperature (RTL12C006-WH only). The contoured seat fits most standard and elongated toilets and includes a lid for added hygiene and aesthetic.

Drive PreserveTech™ products are manufactured with an antimicrobial\* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.

\*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

**EASY-TO-USE CONTROLS**  
Activate and adjust water pressure and select cleansing modes



**HYGIENIC & GENTLE**  
Provides front or back cleansing depending on needs

**EXTRA HEIGHT**  
Contoured seat adds 4" of height to standard or elongated toilets for safer transfer

**FITS MOST STANDARD & ELONGATED TOILETS**



**USER-FRIENDLY**  
Easily connects to water source\*, no electric connection or batteries required



ITEM #	PRODUCT DESCRIPTION	UPC	UOM
RTL12C005-WH	PreserveTech™ Raised Toilet Seat with Bidet, Ambient Water Temperature	822383015941/ea / 50822383015946/cv	1/ea, 1/cv
RTL12C006-WH	PreserveTech™ Raised Toilet Seat with Bidet, Ambient & Warm Water Temperature	822383015958/ea / 50822383015953/cv	1/ea, 1/cv

**PRODUCT SPECS**

Product Dimensions	17.14" (W) x 16.13" (L) x 7.58" (H)
Product Weight	3.5 lb. (RTL12C005-WH)   3.54 lb. (RTL12C006-WH)
Color	White
Weight Capacity	400 lb.
Warranty	One Year Limited**

**BOX SPECS**

Retail Box Dimensions	8.07" (W) x 17.52" (L) x 16.54" (H)
Retail Box Weight	4.5 lb. (RTL12C005-WH)   4.54 lb. (RTL12C006-WH)

\*Initial temperature of water release is contingent upon current temperature of surroundings and will warm based upon your current water setting and supply.

\*\*Refer to instruction guide or website for warranty details.



SCAN FOR BILT. 3D INSTRUCTIONS



POWERED BY BILT INTELLIGENT INSTRUCTIONS\*

FOR USE WITH  
ITEM #RTL12C005-WH

SCAN FOR BILT. 3D INSTRUCTIONS



POWERED BY BILT INTELLIGENT INSTRUCTIONS\*

FOR USE WITH  
ITEM #RTL12C006-WH

