

791 Nitro



- 1. Hand Grip (ea)- 791-01
- 2. Knee Post Adj Pin- 791-02
- 3. Seat Post Release Lever- 791-03
- 4. Knee Pad w/Hardware- 791-04
- 5. Knee Pad Post- 791-05
- 6. Knee Pad Hardware- 791-06
- 7. Front Wheel Hub Cap- 791-07
- 8. Axle Bolt- 791-08
- 9. Right Front Wheel- 791-09
Left Front Wheel- 791-10
- 10. Front Wheel Bearing (pr)- 791-11
- 11. Caster frame Knob- 791-12
- 12. Handlebar- 791-13
- 13. Tiller Adj. Knob- 95021S79907
- 14. Blue Caster Frame Ext.- 791-15BL
Red Caster Frame Ext.- 791-15RD
- 15. Folding Release Lever Assy.- 791-16
- 16. Column to Fork Locking Bolt- 791-17

- 17. Folding Hinge- 791-18
- 18. Column- 791-19
- 19. Fork Blue- 791-20BL
Fork Red- 791-20RD
- 20. Fork Stem Bearing Set- 791-21
- 21. Left Rear Wheel w/Hdwr- 791-22
Right Rear Wheel w/Hdwr- 791-32
- 22. Wheel Frame Bearings (pr)- 791-23
- 23. Rear Axle Assy- 791-24
- 24. Brake Assy- 791-25
- 25. Brake Cable- 791-26
- 26. Left Brake Handle- 791-28
Right Brake Handle- 791-27
- 27. Cable Retaining Clip- 791-29
- 28. Rear Hub Cap- 791-30
- 29. Basket- 791-33