

## Cable Connector for Quattro & ComboCare, 2/PK

HCPCS:

- The Cable Connector connects the lead wire to your electrotherapy device.
- It is designed for use with the Quattro 2, Quattro 2.5 and ComboCare devices manufactured after January 2013.
- Lead wire connects to socket jack in connector, and cable plugs directly into machine
- Cable features 2-pin plug for lead wire pin connector
- For use with lead wire WQ6005PLUG
- For use with Quattro 2, Quattro 2.5 and ComboCare electrotherapy devices manufactured after January 2013
- Quantity: 2 cables/pack

<b>PRODUCTS</b>	<b>Description</b>
DQ8450CPLUG	CABLE CONN FOR DQ8450 & DQ7844 1PKG/