

5' Lead Wire, 2/PK

HCPCS:

- These Replacement Lead Wires easily connect electrodes to the TENS device for effective electrode therapy.
- They are 5 feet long, longer than standard lead wires, and are compatible with most TENS units.
- Connect electrodes to TENS device with 2mm pin connectors for high quality, effective electrode therapy
- Longer than standard lead wires at 5' (63 inches)
- Compatible with most electrotherapy units; NOT compatible with EMPI devices
- Lead wires should be replaced every 4-6 months
- Quantity: 2/package

PRODUCTS	Description
WW5000	5 Foot Lead Wire(2 per pkg) 1EA/