

Model 646 Nebulizer

HCPCS:



- The Model 646 Nebulizer converts liquid into aerosol via high velocity, compressed air stream Used in conjunction with a DeVilbiss brand compressor, it delivers consistent aerosol particles to the lungs of COPD patients
- Used with a DeVilbiss brand compressor for treatment of asthmatic patients
- Designed for comfortable administration of fine mists of uniform density to reach deep into lung capacity
- Unique "T" piece allows patient to leave nebulizer in mouth while exhaling
- Attaches to aerosol face mask or face tent and can be hand operated or compressed air driven with a DeVilbiss brand compressor
- Break-resistant nebulizer that easily disassembles for sterilization by using one of the following methods: autoclaving, boiling or cold sterilization
- Polycarbonate plastic nebulizer may be used with or without mouthpiece (No 25-31-605)
- Comes complete with tubing, activator valve, mouth piece, and extra jet

PRODUCTS

Description

646	Nebulizer
-----	-----------

SPECIFICATIONS

Description

MMAD	5 microns
------	-----------

ACCESSORIES

Description

642-7-605	FINGER VALVE 5 5PK SEPARATE SALE PACKAGE, PKG
-----------	---

25-31-605	MOUTHPIECE W/OUT AIR VALVE 5 PK SS, PKG
-----------	---