

Molded Foam General Use Wheelchair Cushion

HCPCS: E2601-1,E2601-1





FIGURE A

- A unique construction of high-resilient molded foam provides optimal support, stabilization and pressure redistribution
- Designed for low risk patients seeking comfort and support
- Removable, zippered, low-shear stretch cover is vapor permeable and water-resistant
- Cover has a "Non-Slip" vinyl bottom for comfort and safety
- Convenient 12" tie straps secure the cushion to wheelchair
- Meets flammability standard TB117

Description	
W/C Cushion Molded,16x16x2, 1/ea, 6/cv	
W/C Cushion Molded, 20x16x2, 1/ea, 6/cv	
W/C Cushion Molded, 18x16x2, 1/cv, 6/ct	
W/C Cushion Molded, 18x18x2, 1/ea, 6/cv	
W/C Cushion Molded, 20x18, 1/ea, 6/cv	
Description	
18 Months Limited	
250 lbs 300 lbs	
Fluid-Resistant Stretch Nylon	
	W/C Cushion Molded, 16x16x2, 1/ea, 6/cv W/C Cushion Molded, 20x16x2, 1/ea, 6/cv W/C Cushion Molded, 18x16x2, 1/cv, 6/ct W/C Cushion Molded, 18x18x2, 1/ea, 6/cv W/C Cushion Molded, 20x18, 1/ea, 6/cv Description 18 Months Limited 250 lbs 300 lbs