

Titanium Skin Protection and Positioning Gel/Foam Wheelchair Cushion

HCPCS: E2607




FIGURE A



FIGURE B

- The Titanium cushion is specifically designed to offer superior pressure redistribution and positioning
- High-resilient, high-density foam combined with a dual chamber gel bladder encourages proper patient posture through immersion
- Middle contour acts as a built-in leg abductor
- Low-shear, fluid-resistant stretch cover helps to maintain sensitive skin tissue
- Non-skid bottom offers additional stability
- Meets flammability standard TB117

PRODUCTS	Description
FPT-1	TitaniumViscoelasticfoam/GelCush16x16x3
FPT-1818	Titanium Gel/Foam W/C cushion 18X18X3.5
FPT-2	TitaniumViscoelasticfoam/GelCush18x16x3
FPT-2018	Titanium Gel/Foam W/C cushion 20X18X3.5
FPT-4	TitaniumViscoelasticfoam/GelCush22x18x3

SPECIFICATIONS	Description
Warranty	1 Year Limited
Weight Capacity	Up to 20" 275 lbs 22" and Up 500 lbs
Cover Material	Fluid Resistant Stretch Top Non-Skid Bottom
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]