

## Trunk Support

HCPCS: E1399-1, E1399-1



FIGURE A

- Secures the torso at desired angle
- Height, depth, width and angle adjustable
- Padded for cushioned comfort
- Stays open for easy transfers
- Can be mounted towards the inside or outside of handlebar
- Adjustable straps provide a range of support from snug to loose depending on user's capabilities
- In its lowest position, can be used as pelvic positioner
- Adjustments can be made while user is in the gait trainer
- Small trunk support should not be used on TK 3000

PRODUCTS	Description
TK 1080 L	Trunk Support for Moxie, Lg,1/bx
TK 1080 M	Trunk Support for Moxie, Med,1/bx
TK 1080 S	Trunk Support for Moxie, Sm,1/bx