

HCPCS:

PRODUCTS	Description
14001-L	L-Shape Connector, LAL,1/ea
14001T	Tubing for 14001/14002,1/ea
14001-T	T-Connector for LAL, 1/ea
14007-3PORT	Tubing Adapter 3-Port,14029DP,1/ea
14007OR-3	O Ring for 14007, 3/pkg
14013	CPRPlugForLowAirLossMatt,3Port,1 ea
14027-BASE	Replacement Base for 14027, 1/ea
14027-CELL	Repl Bladder - 14027, 1/ea
14028T	Tubing,7", LAL, 1/ea
14028V	OneWayValveForLowAirLossMatt,1/ea
14029-BASE	Base Only for 14029, 1/ea
14029C	Cover for 14029,1/ea
14029-CELL	Bladder Only,14029,1/ea
14029-CELL-1	Cell Only wLazer Holes,14029,1/ea
14029-FFTC	Connector for 14029 Cell, Female, 1/ea
14029-FMTC	Connector for 14029 Cell, Male, 1/ea
14029MN	Mattress Only,14029,3Tube System,1/ea
14030XP-S	Pump Socket 3 Hose Connection,1/ea