

STRENGTHBRACE Universal Knee Compression Sleeve

HCPCS:

- The STRENGTHBRACE Universal Knee Compression Sleeve is the perfect addition for patients who wear a hinged knee brace for compression purposes.
- With a neoprene lining, it's comfortable and prevents the knee brace from sliding or moving.
- Wrap-around style is universal, fitting a wide range of patient sizes
- Neoprene lining keeps knee brace in place
- Decreased material bulk for added patient comfort
- Can be worn on either knee
- To give patients added comfort and support to their knee brace, purchase the STRENGTHBRACE Universal Knee Compression Sleeve.

PRODUCTS

Description

BK2397-2	StrengthBrace Suspension Slv 25EA/CV
----------	--------------------------------------