

StrengthBrace Universal Wrist Brace (Right)

HCPCS:

- Constructed of lightweight, breathable neoprene material, the StrengthBrace Universal Wrist Brace provides support and stabilization to those with wrist pain.
- Designed for those with carpal tunnel syndrome, sprains, strains, tendonitis, wrist pain.
- Universal sizing
- Adjustable dorsal stay
- Lightweight, breathable material
- Easy wrap design for quick fitting
- Inset, stretchable, soft sock for ideal comfort
- Not made with natural rubber latex
- One-year warranty
- PDAC Assigned Code: L3908*

PRODUCTS

Description

BW3908R

STRENGTHBRACE UNIVERSAL WRIST 25EA/CV