

Meridian Alternating Pressure Pad with Electric Pump

HCPCS:



FIGURE A

- The Meridian SatinAir Adjustable Alternating Pressure Pump and Pad System is designed to be placed over an existing mattress and is both easy to use and maintain.
- The compact, quiet pump alternately inflates and deflates air cells in the pressure pad and comes standard with a variable pressure setting for maximum patient comfort and compliance.
- The pump features an output of five liters per minute with a five minute cycle time, and can be mounted on the bed with built-in brackets.
- FEATURES:
 - 5 liter-per-minute pump alternately inflates and deflates air cells over a 5 minute cycle time
 - Variable pressure setting for added patient comfort
 - Pump can be mounted on the bed with built-in bracket
 - Pad includes flaps that wrap around mattress for better stability
 - Does not contain latex
 - Support patients up to 300 lbs.
 - One-year limited warranty on pump, 90-days limited warranty on pad

PRODUCTS

Description

MER-APP2

SATINAIR APP SYSTEM DLX PAD 6EA/CV