

Replacement Cover for APM-8000

HCPCS:

- This is the replacement nylon top cover for the Genesis III 8" alternating pressure mattress.
- It is vapor permeable, low shear, and fluid resistant, protecting the mattress and keeping the patient comfortable.
- Standard replacement nylon top cover protects APM and keeps patients comfortable
- Vapor permeable, low shear, and fluid resistant
- Cover only; does not include mattress or pump
- Designed for use with APM-8CS Pump and Genesis III 8" APM (APM-8000-GBN); may fit other 80" x 36" x 8" mattresses
- Sold individually

PRODUCTS	Description
90179	Replacement Cover for APM-8000 1EA/