

Suction Kit for Heavy-Duty Aspirator (ROS-COMP)

HCPCS:

- This is the Suction Kit for the Roscoe Heavy-Duty Aspirator (ROS-COMP).
- Users can replace all necessary parts in their aspirator conveniently, as it includes a replacement elbow, bacteria filter, reusable collection jar, 6' tubing and two 7" tubing.
- Each item in the kit can also be purchased separately if desired:
 - Elbow - #ROS-610E
 - Bacteria Filter - #ILB-ROS
 - Collection Jar - #1160B-1
 - 6' Tubing - #55N56A
 - 7" Tubing - #TUB-S516
- To keep your heavy-duty aspirator performing optimally, replace its parts regularly and easily with the Roscoe Suction Kit.

PRODUCTS	Description
90386	ROSCOE SUCTION KIT FOR ASP-ROS 1PKG/