

Hydrophobic Filter Tubing Kit for Invacare IRC1135

HCPCS:

- This replacement Hydrophobic Filter Tubing Kit provides everything needed to keep the Invacare IRC1135 Suction Aspirator performing optimally in one purchase.
- It includes a bacteria filter, 16" tubing, and an elbow connector.
- Replacement bacteria filter prevents fluids and other contaminants from entering and damaging suction pump
- Durable and flexible 16" tubing with elbow connector
- Compatible with Invacare IRC1135
- To conveniently replace what you need for your Invacare suction aspirator, purchase the replacement Hydrophobic Filter Tubing Kit.

PRODUCTS

Description

AG402SET

AG IND HYDROPHOBIC FILTER 1EA/