

## Cabinet Filter for Respironics Millennium

HCPCS:

- This replacement cabinet filter is designed for use with Respironics Millennium oxygen concentrator units.
- It processes and filters the air as it passes through the concentrator to provide patient with clean oxygen.
- Cleans bacteria from air as it passes through the oxygen concentrator, providing patient with clean oxygen
- Should be removed and cleaned at least once a week; replace as needed
- Compatible with Respironics Millennium units
- Quantity: 20/pack

PRODUCTS	Description
CF-1030	Cabinet Filter, Millennium 20EA/PKG