

Disposable CPAP Filters for Resironics DreamStation

HCPCS:

- These ultra-fine disposable CPAP filters are essential for your device, as they filter the CPAP airflow of airborne particles and pollen.
- They are designed for use with the Philips Resironics DreamStation Series of CPAP machines.
- Effectively filter airborne particles, dirt, and pollen from your CPAP airflow to ensure clean air
- Ultra-fine design fits easily in device
- Disposable; should be replaced every 2-4 weeks for best device performance
- Fit CPAP devices in the Philips Resironics DreamStation Series
- Dimensions: 2 3/8" L x 1" W
- Quantity: 2 filters per pack
- HCPCS Code: A7038
- For cleaner and healthier CPAP therapy with your Resironics DreamStation device, purchase your Disposable CPAP Filters.

PRODUCTS

Description

CPF-DSD2

DREAMSTATION DISPOSABLE FILTER 2EA/PKG