

Suction Bacteria Filter with 1/4-3/8" Stepped Barb x Barb

HCPCS:



- The replacement Bacteria Filter prevents fluids and other contaminants from entering and damaging your suction pump.
- It is designed for use with Roscoe Suction Kits 50004 and 50006, SUC-COMP, and 601 and 605 Sunrise 7305D, but it can also be compatible with other suction machines.
- Filter connections include 1/4" to 3/8" stepped barbs on both ends
- For use with Roscoe Suction Kits 50004 and 50006, SUC-COMP, and 601 and 605 Sunrise 7305D, but can also be compatible with other suction machines
- Sold individually
- To best protect your suction machine from contamination, purchase your replacement Bacteria Filter.

PRODUCTS

Description

ILB-38B	Suction Filter,1/4-3/8 1EA/
---------	-----------------------------