

Replacement Hand Pendant for

HCPCS:

- This is the replacement Hand Pendant for the ProBasics 42" Bariatric Bed (#PB42BARBED).
- It is only the pendant and does not include motor or any other parts of the bed.
- Replacement pendant for 42" Bariatric Bed (#PB42BARBED)
- Allows patient or caregiver to easily adjust head and foot sections as well as bed height
- Sold individually
- Purchase the replacement Hand Pendant for the 42" Bariatric Bed to continue easily adjusting patient positions and bed height.

PRODUCTS

Description

PBHB4PNDC

Replacement Hand Pendant for 1EA/