

Piston Cup for Thomas 2619, 2621, and 2639 Compressors

HCPCS:

- This is the replacement piston cup designed for use with Thomas 2619, 2621 and 2639 Compressors.
- The piston cup seal should be replaced regularly for optimal air compressor performance.
- Essential for the operation of compressor pump
- Fits Thomas 2619, 2621 and 2639 Compressors
- To keep your Thomas 2619, 2621 or 2639 Compressor performing optimally, purchase the replacement piston cup.

PRODUCTS	Description
PC-19	Piston Cup 1EA/