

Bestcare Folding Electric Patient Lift

HCPCS:

- The BestLift PL400EF (Drop Ship)
- This versatile lift represents a new standard in full body lifts for individuals who travel or have storage issues in their home.
- It is a unique lifter with a small footprint and foldable design that assembles or folds in 3 simple steps allowing for easy storage or transport.
- As with all Bestcare lifts, this new Bestcare lifter continues to insure maximum safety and comfort.
- Main features:
- Foldable / Portable / Transportable
- Smaller and lighter
- Fits into most car trunks
- 400lb capacity
- 6-point spreader bar (2 point optional)
- New Bestcare Control system

PRODUCTS	Description
PL400EF	BESTCARE ELECTRIC FOLDING PATIENT LIFT