

Carex OTC Finger Pulse Oximeter

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- The Carex Finger Pulse Oximeter allows you to monitor and maintain pulse and oxygen saturation levels for optimal health.
- Accurate, fast oxygen level and heart rate measuring technology.
- Ideal for sport enthusiasts! No prescription required.
- Approved for over-the-counter retail sale
- Large easy-to-read display
- Suitable for pediatric and adult use
- Soft silicone padding for a comfortable fit
- One-button operation
- Up to 30 hours of use
- Auto-off function
- Low battery indicator

PRODUCTS	Description
FGR41800 0000	CAREX FINGERTIP PULSE OXIMETER 100EA/CV