

## Salter PRO2 Check Elite Oxygen Concentration Indicator

HCPCS:

- The Salter PRO2 Check Elite Oxygen Concentration Indicator accurately indicates the O2 percentage, output pressure, and output flow from concentrators.
- It makes it easy to check the quality of oxygen, connecting to concentrators with standard tubing and presenting results clearly on LCD display.
- Includes three operation modes: pressure (measured in PSIG or kPa), flow (measured in LPM), and oxygen concentration (measured in percent)
- Calibrated out of box; calibration is not affected by temperature or altitude to ensure accuracy
- Connects to concentrator with standard oxygen supply tubing
- Easy-to-read LCD display
- Operates on 9 volt battery
- Automatically shuts off after 5 minutes to preserve battery
- Compact and portable; comes with carrying case

### PRODUCTS

### Description

---

PRO2 ELITE-1	SALTER PRO2 CHECK ELITE 1EA/
--------------	------------------------------

---