

PROBASICS REPLACEMENT SEAL

HCPCS:

- Roscoe Zzz-Mask Full Face Replacement Seal, Medium
- About The Mask
- The Roscoe Zzz-Mask Full Face CPAP Mask with Comfort Cushion is made of high-quality clear silicone with a virtually indestructible frame.
- The lightweight mask comes complete with headgear and dual 360° swivel port which allows for free movement and easy connection at the tubing attachment point.
- The mask is dishwasher safe and comes complete with headgear
- Features
- Made of high-quality, lightweight, clear silicone with a flexible design
- Dual 360° swivel port at the tubing attachment point
- Air vent holes and one-way valve design help to reduce CO2 retention
- Made with half the parts of popular brands, meaning less patient confusion
- Headgear with hook-and-loop fasteners provide maximum versatility and movement

PRODUCTS	Description
PB780M	PROBASICS REPLACEMENT SEAL 100EA/CV