

## PB Style Chin Strap

HCPCS:

- The PB Style Chin Strap prevents air loss during CPAP therapy with a nasal mask by ensuring the patient's mouth stays closed throughout the night.
- In turn, it keeps patient comfortable by reducing side effects such as dry mouth and air swallowing.
- Comfortably keeps mouth shut by supporting the patient's chin
- Easily adjustable with a single hook-and-loop strap that fastens at the top of the head
- Comparable to Puritan Bennett #133318
- Ruby red
- Sold individually
- For a comfortable way to ensure continual oxygen therapy for patients throughout the night, purchase the PB Style Chin Strap.

<b>PRODUCTS</b>	<b>Description</b>
ROS-133318	PB Style Chin Strap 1EA/