

## 3 Point Chin Strap, Large

HCPCS:

- This ruby-style chin strap is a simple and comfortable solution for those who don't want to use a complicated chin strap for their CPAP or BiPAP mask.
- It is made of a breathable material with adjustable straps located at both the crown of the head and at the base of the neck and are made with Velcro material.
- It also aids in keeping your mouth closed during CPAP or BiPAP therapy, and may prevent the patient from having to wear a full face CPAP mask.
- Compare to PureSom Fusion CPAP chinstrap TMS-09L.
- Size: large.
- Sold as each.

<b>PRODUCTS</b>	<b>Description</b>
ROS-T09L	3 Point Chin Strap, Large, 1EA/