

ELECTRIC HOSPITAL BED ENDS 1EA/

HCPCS:

PRODUCTS	Description
15033BV-PKG-MED	UltraLtFullElecPkg, 2/PAL
HBFEBE	Full-Electric Two-Motor Bed 1EA/
HBFESS	Full-Electric Two Motor 1EA/
PB3000	BARIATRIC BED, 42", HEAD SPRNG 1EA/
PB3010	BARIATRIC BED, 42", FOOT SPRNG 1EA/
PB3301	BARIATRIC BED, 42", BED ENDS 1EA/
PBFEB	2-MOTOR FULL-ELEC BED 1EA/
PBFEB-FRARBPKG	2-MOTOR FULL-ELEC BED PACKAGE 1EA/
PBFEBLO-HRARBPKG	2-MOTOR FULL-ELEC BED LOW PKG 1EA/