

SALTER HIGH FLOW CANNULA

HCPCS:

- High Flow O2 Nasal Cannula features include:
- Delivers oxygen flows up to 15 LPM
- Special three channel headset tubing with enhanced reservoir face piece
- Unique, rectangular style head set tubing slide bolo
- Seamless, extra soft, reservoir face piece anatomically contoured to fit upper lip
- Bright green supply tubing from the cannula to oxygen source

PRODUCTS

Description

1600HF-50

SALTER HI FLOW CANNULA W/50' TUB 10EA/CV