

EQ 02 CANNULA 6FT 12EA/CV

HCPCS:

PRODUCTS	Description
1056	SALTER CANNULA, W/SMOOTH BORE 1CV/
1600-14	SALTER CANNULA W/14' TUBE 50EA/CV
1600-25	SALTER CANNULA W/25' TUBE 25EA/CV
1600-4	SALTER CANNULA W/4' TUBE 50EA/CV
1600HF-14-25	SALTER HI FLOW CANNULA W/14' TUB 25EA/CV
1600HF-25	SALTER HI FLOW CANNULA W/25' TUB 10EA/CV
1600HF-4-25	SALTER HI FLOW CANNULA W/4' TUBE 25EA/CV
1600HF-50	SALTER HI FLOW CANNULA W/50' TUB 10EA/CV
1600HF-7-25	SALTER HI FLOW CANNULA W/7' TUBE 25EA/CV
1600TLC	SALTER TLC CANNULA W/7' TUBE 25EA/CV
1615-7 50	SALTER LABS, INTERMEDIATE
1616-4	SALTER MICRO CANNULA W/4' TUBE 25EA/CV
16SOFT-25-25	SALTER SOFT CANNULA W/25' TUBE 25EA/CV
16SOFT-4-50	SALTER SOFT CANNULA W/4' TUBE 50EA/CV
16SOFT-7-50	SALTER SOFT CANNULA W/7' TUBE 50EA/CV
190	Comfort Plus Cannula Harness, 1CV/
195-W	Adult Cannula with 25 ft. 1CV/
197-W	Comfort Plus Cannula, 7 ft. 1CV/
198-W	Adult Cannula with 14 ft. kink 1CV/
27-W	Kink-resistant tubing, 25 ft.,
4804	SALTER DEMAND CANNULA W/4' TUBE 25EA/CV
4807	SALTER DEMAND CANNULA W/7' TUBE 25EA/CV
556	WESTMED CMFT SOFT PLUS CANNULA 1CV/
588	Adult Nasal Cannula with 1EA/
CAN-PEDROS	PEDIATRIC O2 CANNULA W/7' TUBE 50EA/CV
CAN-ROS	OXYGEN SOFT CANNULA W/O TUBE 50EA/CV
CAN-ROS25	OXYGEN SOFT CANNULA W/ 25' TUBE 25EA/CV
CAN-ROS4	OXYGEN SOFT CANNULA W/ 4' TUBE 50EA/CV
CAN-ROS7	OXYGEN SOFT CANNULA W/ 7' TUBE 50EA/CV
CON 300	PREMIUM OXYGEN SWIVEL ADAPTER 100EA/CV
FGR10400 0000	Comfort SoftPlus Cannula 7' tube 24EA/CV
HUM 300	PREMIUM BUBBLE HUMIDIFIER-6 PSI 50EA/CV

PRODUCTS

Description

HUM HF300	PREMIUM HI-FLO BUBBL HUMIDIFIER 25EA/CV
MASK 300A	PREMIUM ADULT AEROSOL MASK 50EA/CV
MASK 300AO	PREMIUM OXYGEN MASK W/ 7' TUBE 50ES/CV
MASK 300P	PREMIUM PEDIATRIC AEROSOL MASK 50EA/CV
NEB KIT 300	PREMIUM NEBULIZER W/T & 6" RES 50EA/CV
NEB KIT 300AM	PREMIUM NEBULIZER W/ADULT MASK 50EA/CV
NEB KIT 300PM	PREMIUM NEBULIZER W/PED MASK 50EA/CV
SOFT 300	PREMSOFT CANNULA W/O TUBE 50EA/CV
SOFT 307P	PREMSOFT PED CANNULA W/ 7' TUBE 50EA/CV
SOFT 325	PREMSOFT CANNULA W/ 25' TUBE 25EA/CV
SOFT HF307	PRSOFT HI-FLO CANNULA W/7' TUBE 25EA/CV
SOFT HF314	PRSOFT HI-FLO CANNULA W/14' TUBE 25EA/CV
SOFT HF325	PRSOFT HI-FLO CANNULA W/25' TUBE 10EA/CV
TUB 307C	PREMIUM 7' OXYGEN TUBING-CLEAR 50EA/CV
TUB 325C	PREMIUM 25' OXYGEN TUBING-CLEAR 25EA/CV
TUB 325G	PREMIUM 25' OXYGEN TUBING-GREEN 25EA/CV
TUB 350C	PREMIUM 50' OXYGEN TUBING-CLEAR 20EA/CV
TUB 350G	PREMIUM 50' OXYGEN TUBING-GREEN 20EA/CV
WTR 300	PREMIUM IN-LINE WATER TRAP 25EA/CV