

GEN PURPOSE EQUIP BAGS 1.5 MIL 30X42

HCPCS:

PRODUCTS	Description
BAG-211330R	CONCENTRATOR EQUIP BAGS 1 MIL 21X13X30
BAG-2430R	GEN PURPOSE EQUIP BAGS 1.5 MIL 24X30
BAG-2636R	GEN PURPOSE EQUIP BAGS 1.5 MIL 26X36
BAG-3042R	GEN PURPOSE EQUIP BAGS 1.5 MIL 30X42